



Defender BELMED INC.

Automatic Changeover Manifold



Defender Automatic Changeover Manifold



The Belmed Defender Automatic Changeover Manifold is the perfect manifold solution for your dental nitrous oxide delivery system. The Defender can accommodate two cylinders for both oxygen and nitrous oxide. This manifold is fully automatic and will switch from "In- Use" to "Reserve" tank without fluctuations in the gas delivery. The "Empty" indicator will remain lit until a new cylinder is connected and the appropriate pressure is detected by the manifold. The system is auto-sensing...no need to manually reset the manifold after a cylinder change.

- Easy to read LED display clearly shows gas pressures and cylinder status
- Automatic detection of cylinder pressure, no manual reset is required
- Medical grade sensors that give you the reliability you need
- Intergrated Dry Contacts for additional remote monitoring
- Complies with NFPA 99

Packages



Packages with Pre-Install Kit

Part Number	Description	MSRP
A100	Defender Automatic Changeover Manifold with Wall Alarm & Pre-Install Kit – Includes: Automatic Changeover Manifold, Wall Alarm, Regulators, Hoses, Tank Restraints, and Pre-Install Kit	\$4,900.00
A101	Defender Automatic Changeover Manifold with Desk Alarm & Pre-Install Kit – Includes: Automatic Changeover Manifold, Desk Alarm, Regulators, Hoses, Tank Restraints, and Pre-Install Kit	\$4,900.00

Packages without Pre-Install Kit

Part Number	Description	MSRP
A120	Defender Automatic Changeover Manifold with Wall Alarm & Pre-Install Kit – Includes: Automatic Changeover Manifold, Wall Alarm, Regulators, Hoses, and Tank Restraints	\$4,680.00
A121	Defender Automatic Changeover Manifold with Desk Alarm & Pre-Install Kit – Includes: Automatic Changeover Manifold, Desk Alarm, Regulators, Hoses, and Tank Restraints	\$4,680.00

Pre-Install Kit

Part Number	Description	MSRP
A130	Defender Pre-Install Kit – Includes: Pipeline Connectors, Cable, and Gas Identification Labels.	\$220.00

Zone Valve



Part Number	Description	MSRP
4054	Zone Valve – Includes: Zone Valve for Defender Automatic Changeover Manifold and Cable	\$1,350.00

Additional Accessories

Part Number	Description	MSRP
4024-0050	Defender O2 Pipeline Connector	\$72.00
4024-0051	Defender N2O Pipeline Connector	\$72.00



OXYGEN

NITROUS OXIDE


Part Number	Description	MSRP
4022BA-28	O2 Gas Identification Label (6 Pkg)	\$10.00
4022BA-30	N2O Gas Identification Label (6 Pkg)	\$10.00

Part Number	Description	MSRP
4024-0002-0001	Defender Transformer	\$80.00



Part Number	Description	MSRP
4024-0060-0001	Defender Cable (100ft)	\$80.00
4054-0002-0001	Defender Zone Valve Cable	\$2.50/ft

Part Number	Description	MSRP
8138-N	Tank Restraint	\$112.00



Part Number	Description	MSRP
8005X4	O2 DISS Nut x DISS Nut Hose Assembly (4 Feet)	\$86.00
8006X4	N2O DISS Nut X DISS Nut Hose Assembly (4 Feet)	\$86.00

Part Number	Description	MSRP
6001	O2 Single High Pressure Gauge – 50 PSI Pre-Set	\$326.00
6002	N2O Single High Pressure Gauge – 50 PSI Pre-Set	\$326.00



Order:
orders@belmedinc.com
Fax: 717-246-7586

Mail To:
PO Box 247
Red Lion, PA 17356

Contact Us:
Local: 717-246-5500
Toll Free: 888-723-5893
bmi@belmedinc.com

Ship To:
887 Delta Rd
Red Lion, PA 17356

Tech Support:
help@belmedinc.com