

HEADREST SELECTIONS, ALL MODELS




Lever Lok


Dual Articulating Headrests
Release lever, position, and lock
into place. Allows headrest to
be positioned over the back for
smaller patients and articulates
over 90 degrees for airway
management. Available with and
without forearm supports.



LL-69 Lever
Pillow Prosthetic
P/N 67-1569 Putty Trim
P/N 67-1569B Black Trim



LL-62 Lever
Pillow Prosthetic
w/Forearm Supports
P/N 67-1562 Putty Trim
P/N 67-1562B Black Trim



LL-63 Lever
Flat w/Magnetic Pillow
P/N 67-1563 Putty Trim
P/N 67-1563B Black Trim



Pinch Lok

Dual Articulating Headrests
Pinch tabs for positioning and release
to lock. Available with and without
forearm supports.



LL-8 Lever
Tapered Rectangle
P/N 67-1508 Putty Trim
P/N 67-1508B Black Trim



LL-10 Lever
Omega Shaped
P/N 67-1560 Putty Trim
P/N 67-1560B Black Trim




LL-64 Lever
Flat w/Magnetic Pillow &
Forearm Supports
P/N 67-1564 Putty Trim
P/N 67-1564B Black Trim



PL-49 Pinch LOK
Pillow Prosthetic
P/N 67-1549 Putty Trim
P/N 67-1549B Black Trim



PL-42 Pinch LOK
Pillow Prosthetic w/
Forearm Supports
P/N 67-1542 Putty Trim
P/N 67-1542B Black Trim



PL-53 Pinch LOK
Flat w/Magnetic Pillow
P/N 67-1553 Putty Trim
P/N 67-1553B Black Trim



PL-54 Pinch LOK
Flat w/Magnetic Pillow &
Forearm Supports
P/N 67-1554 Putty Trim
P/N 67-1554B Black Trim