













## Whicam Story3 CS – Wireless Intraoral Camera

# Enjoy Wireless & Auto Focus

## Wireless Autofocus Intraoral Camera with Mouse

### 1 Features of Whicam Story3 CS

-  **High speed digital wireless**  
High speed digital wireless provide a clearer live image without interference.
-  **Auto focus**  
Real Autofocus function, semi-focus also available.
-  **HD image**  
High definition, 1280x720 Image
-  **Gyro mouse built-in**  
Useable as a mouse.
-  **Adjustable head angle**  
Three adjustable angles to reduce arm fatigue.
-  **280° twist head**  
Easy to maneuver in patients mouth.
-  **Larger battery capacity**  
User replaceable and longer battery life, adding more hours of usage.
-  **Software compatibility**  
Integrates with almost all softwares or we will provide a free imaging software.
-  **OLED screen**  
Bright and clear OLED screen indicates various device status.
-  **Adjustable LED brightness**  
LED intensity adjustment button, you can adjust brightness to three intensity levels.



### 2 Wireless Receiver

- USB power operation  
It doesn't need other electric power supply.  
The receiver works with USB electric power only.
- Compact design  
Compact receiver saves space.
- Unlimited receiver channels  
Digital wireless type supports unlimited channels.  
User can use many receivers with only one camera.

### 4 Minimum System Requirements

|            |                                     |
|------------|-------------------------------------|
| PC system: | 2.0 Ghz CPU                         |
| OS:        | Windows 7 or 8 Windows XP Not HIPAA |
| USB port:  | USB 2.0                             |
| RAM:       | 4 GB RAM                            |

### 3 Specifications

|                           |                   |                                 |
|---------------------------|-------------------|---------------------------------|
| Camera<br>(included head) | Operation voltage | 3.7V                            |
|                           | Measurement       | Approx. 266 x 31 x 37 mm        |
|                           | Weight            | Approx. 140g                    |
|                           | Resolution        | HD 720p                         |
|                           | Sensor/Lens type  | CMOS / Liquid lens              |
| Charger                   | Measurement       | Approx. 169 x 148 x 89 mm       |
|                           | Weight            | Approx. 170 g                   |
|                           | Operation voltage | 5V                              |
| Disposable cover          | Dimensions        | Approx. 32 x 100 mm             |
|                           | Quantity          | Approx. 100 pcs by weight x 2ea |
|                           | Raw material      | Polyethylene                    |
| Adapter                   | Input voltage     | AC 100 ~ 240 V                  |
|                           | Output voltage    | DC 5V 2A                        |
|                           | UL number         | E333054                         |

### 5 Integration

Integrates with almost ALL software's; or we provide a free Imaging software.

