accuc@ld

ARG31PVBIADADL2B



32" x 18.5" x 22.63" (H x W x D)

Performance series Pharma-Vac 2.83 cu.ft. ADA height glass door commercial all-refrigerator for the display and refrigeration of vaccines; designed for freestanding or recessed installation, with factory-installed data logger

Highlights:

32" height fits under lower ADA compliant counters Recessed installation capable for use in 19" wide spaces Intelligent microprocessor digital temperature controller

Product Features:

| Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC guidelines |
|--|
| 32" height is designed to fit under ADA compliant counters |
| The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period |
| Versatile design allows recessed installation under counters or freestanding use |
| Accucold DL2B data logger is NIST calibrated in our ISO/IEC 17025:2017 laboratory, with user-set logging intervals, a large LCD readout of the current and high/low temperature, and virtually unlimited storage through internal memory and a convenient USB drive; for additional information on this logger's features, view the DL2BKit on our website |
| Meets UL-471 commercial standards |
| Intuitive controls let you adjust the interior temperature between +2°C to +8°C |
| Visual and audible alarms alert if the interior temperature rises or falls out of set range |
| Electronic controller is equipped with a volt-free contact for alarm forwarding to an external remote system (contacts located in the external rear of the unit) |
| |

accuc@ld

ARG31PVBIADADL2B Specifications:

| Overview | |
|---|--|
| Height of Cabinet | 32.0" (81 cm) |
| Width | 18.5" (47 cm) |
| Width with Door Open | 19.0" (48 cm) |
| Depth | 22.63" (57 cm) |
| Depth with Handle | 24.38" (62 cm) |
| Depth with door at 90° | 40.0" (102 cm) |
| Capacity | 2.83 cu.ft. (80 L) |
| Defrost Type | Automatic |
| Door | Glass |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 0.75 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 102.0 lbs. (46 kg) |
| Shipping Weight | 113.0 lbs. (51 kg) |
| Parts & Labor Warranty | 2 Years |
| 0 14/ / | |
| Compressor Warranty | 5 Years |
| Refrigerator Features | 5 Years |
| | 5 Years |
| Refrigerator Features | |
| Refrigerator Features Door Swing | RHD |
| Refrigerator Features Door Swing Reversible | RHD Factory Reversible |
| Refrigerator Features Door Swing Reversible Shelf Type | RHD Factory Reversible Wire |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty | RHD Factory Reversible Wire 4 |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves | RHD Factory Reversible Wire 4 Yes |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type | RHD Factory Reversible Wire 4 Yes Digital |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type | RHD Factory Reversible Wire 4 Yes Digital Interior |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type | RHD Factory Reversible Wire 4 Yes Digital Interior R600a |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI Low Side PSI | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI Low Side PSI Interior Light | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 90.0 Yes |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI Low Side PSI Interior Light Temperature Range | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 90.0 Yes |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI Low Side PSI Interior Light Temperature Range Dimensions | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 90.0 Yes 2°C to 8°C |
| Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount High Side PSI Low Side PSI Interior Light Temperature Range Dimensions Interior Height | RHD Factory Reversible Wire 4 Yes Digital Interior R600a 1.41oz. 280.0 90.0 Yes 2°C to 8°C |

| Compressor Step Width | 15.25" (39 cm) |
|-----------------------|----------------|
| Compressor Step Depth | 5.5" (14 cm) |



