

EL11LT



33.5" x 28.5" x 27.25" (H x W x D)

Commercial -45° C capable chest freezer with digital thermostat and 4.8 cu.ft. capacity

Highlights:

Capable of -45° C operation for safe storage of vaccines and other temperature sensitive materials

UL-S listed to NSF-7 commercial standards

Digital thermostat with external controls allows easy and precise temperature monitoring

Product Features:

-45° capable	Low temperature operation is ideal for select vaccines, pharmaceuticals, and frozen goods
Commercially approved	UL-S listed to NSF/ANSI-7 standards
Extra thick insulation	4 1/2" thick walls for added efficiency and stable cooling
Factory installed lock	Keyed lock on the handle for reliable security
Counterbalanced hinges	Easy and safe loading and unloading with a lid that remains open when lifted past 45° angles
Casters	Four solid wheels included for optional mobility
Digital thermostat	Electronic controls for precise temperature management, with temperature readout and controls externally located for added convenience
Manual defrost	Static cooling for a stable interior
100% CFC Free	Environmentally friendly design with no ozone-damaging chemicals

EL11LT Specifications:

Overview	
Height of Cabinet	33.5" (85 cm)
Width	28.5" (72 cm)
Depth	27.25" (69 cm)
Depth with Handle	28.75" (73 cm)
Capacity	4.8 cu.ft. (136 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	7.2
Voltage/Frequency	115 V AC/60 Hz
Weight	110.0 lbs. (50 kg)
Shipping Weight	115.0 lbs. (52 kg)
Estimated Time to Ship	Next day shipping
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Freezer	
Reversible	No
Thermostat Type	Digital
Refrigerant Type	R290a
Refrigerant Amount	1.94 oz.
High Side PSI	342.0
Low Side PSI	174.0

