

# EL21LT



33.5" x 41.5" x 27.25" (H x W x D)

Commercial -45° C capable chest freezer with digital thermostat and 8.1 cu.ft. capacity

### Highlights:

Capable of -45° C operation for safe storage of vaccines and other temperature sensitive materials

UL-S listed to NSF-7 commercial standards

Digital thermostat with external controls allows easy and precise temperature monitoring

### Product Features:

Commercially approved	UL-S listed to NSF/ANSI-7 standards
-45° capable	Low temperature operation is ideal for select vaccines, pharmaceuticals, and frozen goods
Extra thick insulation	4 1/2" thick walls for added efficiency and stable cooling
Factory installed lock	Keyed lock on the handle for reliable security
Digital thermostat	Electronic controls for precise temperature management, with temperature readout and controls externally located for added convenience
Access port	Factory installed probe hole allows users to install their own monitoring equipment
Casters	Four solid wheels included for optional mobility
Counterbalanced hinges	Easy and safe loading and unloading with a lid that remains open when lifted past 45° angles
Manual defrost	Static cooling for a stable interior

## EL21LT Specifications:

Overview	
Height of Cabinet	33.5" (85 cm)
Width	41.5" (105 cm)
Depth	27.25" (69 cm)
Capacity	8.4 cu.ft. (238 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	7.2
Voltage/Frequency	115 V AC/60 Hz
Weight	130.0 lbs. (59 kg)
Shipping Weight	135.0 lbs. (61 kg)
Estimated Time to Ship	Next day shipping
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Freezer	
Thermostat Type	Digital
Refrigerant Type	R290a
Refrigerant Amount	2.29 oz.
High Side PSI	342.0
Low Side PSI	174.0

