

PHC115G



63.75" x 23.63" x 24.13" (H x W x D)

Single chamber warming cabinet with glass door, stainless steel construction, digital thermostat and lock

Highlights:

Stainless steel construction inside and out for durable use and easy cleaning

Digital controls conveniently located for easy adjustment and temperature viewing

Full length glass door keeps complete interior in view for easy inventory check

Product Features:

Digital thermostat	Electronic controls are conveniently located for easy temperature adjustment	
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature	
120 to 150°F Temperature Range	Designed and set for hospital blankets	
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior	
Professional handle	Stainless steel handle completes the professional look	
Cantilevered shelving system	Easily adjustable to accommodate virtually any size	
Stainless steel interior	Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning	
Fully finished stainless steel cabinet	Complete stainless steel construction allows freestanding use	
Factory installed lock	Keyed lock for a secure interior	





PHC115G Specifications:

Overview		
Height of Cabinet	63.75" (162 cm)	
Width	23.63" (60 cm)	
Depth	24.13" (61 cm)	
Depth with Handle	25.25" (64 cm)	
Depth with door at 90°	46.25" (117 cm)	
Capacity	11.0 cu.ft. (311 L)	
Door	Glass	
Cabinet	Stainless Steel	
US Electrical Safety	ETL	
Canadian Electrical Safety	ETL-C	
Amps	5.5	
Voltage/Frequency	115 V AC/60 Hz	
Weight	145.0 lbs. (66 kg)	
Shipping Weight	150.0 lbs. (68 kg)	
Estimated Time to Ship	Next day shipping	
Parts & Labor Warranty	1 Year	
Warming Cabinet		
Door Swing	RHD	
Reversible	No	
Adjustable Shelves	Yes	
Shelf Type	Wire	
Shelf Qty	4	
Thermostat Type	Digital	
Interior Height	50.0" (127 cm)	
	40 5" (50 200)	
Interior Width	19.5" (50 cm)	
Interior Width Interior Depth	19.5" (50 cm) 18.75" (48 cm)	





