

SPRF36LCART



38.5" x 39" x 21.5" (H x W x D)

304 grade stainless steel fully assembled rolling cart RPC with SPRF36 portable 12V/24V vaccine cooler with factory-installed lock

Highlights:

Mobile vaccine station with wheeled stainless steel cart and SPRF36 portable refrigerator-freezer for use indoors or outdoors

SPRF36 refrigerator-freezer can operate as a refrigerator (2-8°C) or freezer (-15°C) for storing medical supplies

Weatherproof 304 grade stainless steel constructed cart for added durability in indoor and outdoor medical environments

Product Features:

Fully assembled design	Cart ships fully assembled for easy use. Just install the handle and place the refrigerator on top
Portable design	Combines our stainless steel wheeled cart with SPRF36 portable refrigerator-freezer for a completely mobile vaccine station
Stainless steel cart	Weatherproof 304 grade stainless steel construction for added durability in indoor and outdoor medical environments
SPRF36 refrigerator-freezer	Portable 12V/24V cooler capable of operating as a refrigerator (2-8°C) or freezer (-15°C) for unparalleled storage flexibility including samples, vaccines, cold packs, and other healthcare supplies
Factory-installed wheels	Four wheels come pre-installed on the cart for easy mobility
Stainless steel cabinet	Non-refrigerated cabinet with fixed shelf and magnetically locking door provides ample storage of medical supplies and other materials
Exterior stainless steel shelf	Side shelf affords extra storage as needed
12V/24V DC hookup	Plug the SPRF36 into the lighter socket of any vehicle or 24V battery
AC Adapter included	We include an AC adapter for users to plug the SPRF36 into a standard 110-120V wall socket for stationary use

SPRF36LCART Specifications:

Overview	
Height of Cabinet	38.5" (98 cm)
Width	39.0" (99 cm)
Depth	21.5" (55 cm)
Capacity	1.0 cu.ft. (28 L)
Defrost Type	Manual
Door	Gray
Cabinet	Stainless Steel
Amps	4.16
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Shipping Weight	150.0 lbs. (68 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Refrigerator-Freezer	
Door Swing	Lift-Up
Refrigerator - Interior Height	15.0" (38 cm)
Refrigerator - Interior Width	13.25" (34 cm)
Refrigerator - Interior Depth	9.75" (25 cm)
Refrigerator - Temperature Range	2-8°C
Refrigerator - Refrigerant Type	R134a
Refrigerator - Refrigerant Amount	1.58 oz.
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	1.58 oz.
High Side PSI	212.05
Low Side PSI	16.68
Condensor Location	Side of Unit
Freezer - Temperature Range	down to -15°C

