

VLT650ADA



32" x 23.63" x 23.5" (H x W x D)

Built-in undercounter ADA compliant laboratory freezer capable of -35° C (-31° F) operation

Highlights:

Commercially listed to NSF-7 standards and meets UL-471

Capable of operating at -35° C

Precise temperature control with a digital thermostat

Product Features:

Commercially approved	UL-S listed to NSF-7 commercial standards and meets UL commercial standard 471
ADA compliant	32" height allows installation under lower ADA compliant counters
-35° C (-31° F) capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Built-in capable	Front-breathing system allows the unit to be used built-in
Factory installed lock	Keyed lock offers convenient security
Three slide-out drawers	Conveniently organize your items while keeping the temperature stable when the door is opened
Digital thermostat	Electronic controls for precise temperature management
Alarm with temperature display	Audible alarm sounds to warn of temperature variations with a digital readout on the thermometer to provide the interior temperature
Manual defrost	Static cooling system with manual defrost operation

VLT650ADA Specifications:

Overview	
Height of Cabinet	32.0" (81 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Capacity	3.5 cu.ft. (99 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Voltage/Frequency	115 V AC/60 Hz
Shipping Weight	135.0 lbs. (61 kg)
Estimated Time to Ship	5 to 7 days
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	2.76 oz.
Interior Drawers	Yes
Drawers Qty	3

