accuc@ld

ARG1PV



ACCUCOLD DIVISION, FELIX STORCH, INC. | ISO 9001:2015 Certified

21.5" x 17.5" x 19.75" (H x W x D)

Performance series Pharma-Vac 1 cu.ft. compact glass door all-refrigerator for vaccine storage

Highlights:

Intelligent microprocessor digital temperature controller

Improved temperature accuracy with a separate temperature probe for continuous product readout to simulate stored product temperature

Additional sensor reads interior air temperature for superior temperature control

Product Features:

Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC guidelines
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Adjustable temperature	Intuitive controls let you adjust the interior temperature between +2 to +8°C
Temperature alarm	Visual and audible alarms alert if the interior temperature rises or falls out of set range
Remote alarm contacts	Electronic controller is equipped with a volt-free contact for alarm forwarding to an external remote system (contacts located in the external rear of the unit)
Automatic defrost with adjustable cycle defrost	Adjustable cycle defrost operation ensures the unit avoid temperature spikes and ice buildup
Access port	Includes a pre-installed 1/2" probe hole in the right interior wall to allow easier connection to user-provided monitoring equipment
Open door alarm	Alarm will sound if the door is left ajar for over 4 minutes (adjustable from 1-99 minutes)
Audible power failure alarm	Battery-powered alarm sounds if the power connection is interrupted



accuc@ld

ARG1PV Specifications:

Overview	
Height of Cabinet	21.5" (55 cm)
Width	17.5" (44 cm)
Width with Door Open	18.0" (46 cm)
Depth	19.75" (50 cm)
Depth with Handle	21.5" (55 cm)
Depth with door at 90°	36.0" (91 cm)
Capacity	1.0 cu.ft. (28 L)
Defrost Type	Automatic
Door	Glass
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.6
Voltage/Frequency	115 V AC/60 Hz
Weight	63.0 lbs. (29 kg)
Shipping Weight	73.0 lbs. (33 kg)
Estimated Time to Ship	1 to 2 days
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
Compressor Warranty Refrigerator Features	5 Years
	5 Years
Refrigerator Features	
Refrigerator Features Door Swing	RHD
Refrigerator Features Door Swing Reversible	RHD Yes
Refrigerator Features Door Swing Reversible Shelf Type	RHD Yes Wire
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty	RHD Yes Wire
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves	RHD Yes Wire 3 Yes
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type	RHD Yes Wire 3 Yes Digital
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type	RHD Yes Wire 3 Yes Digital Interior
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type	RHD Yes Wire 3 Yes Digital Interior R600a
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount	RHD Yes Wire 3 Yes Digital Interior R600a 0.71oz.
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount Interior Light	RHD Yes Wire 3 Yes Digital Interior R600a 0.71oz. Yes
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount Interior Light Temperature Range	RHD Yes Wire 3 Yes Digital Interior R600a 0.71oz. Yes
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount Interior Light Temperature Range Dimensions	RHD Yes Wire 3 Yes Digital Interior R600a 0.71oz. Yes 2 to 8°C
Refrigerator Features Door Swing Reversible Shelf Type Shelf Qty Adjustable Shelves Thermostat Type Fan Type Refrigerant Type Refrigerant Amount Interior Light Temperature Range Dimensions Interior Height	RHD Yes Wire 3 Yes Digital Interior R600a 0.71oz. Yes 2 to 8°C



