

### **Snap and Evac**

Just snap in handpiece of your choice

### **Efficient Workflow**

Eliminates need for OR personnel to hold tubing

# **Transparent Capture Port**

Improves visibility of the surgical site

### **Cost-Effective**

Economical option for surgical smoke evacuation

# **Lightweight Tubing**

QapEvac"

Showing the thoropolar transpace.

#### **Product Ordering Information**

Description	Catalog Number	Quantity
SnapEvac <sup>™</sup> Designed using ABS plastic and a polycarbonate capture port with 10 ft (3 m) tubing; 7/8 in (22 mm) cuff.	SnapEvac20	20/Case

<sup>\*</sup>ESU pencil and handpiece both sold separately. SnapEvac may not fit all ESU pencils and handpieces.



# PenAdapt® Surgical Smoke Evacuation Adapter

## **Universal Design**

Pliable silicone material allows for proper fit with most electrosurgical pencil brands.

# **Angled Back End Clasp**

Reduces tension placed on surgeon's wrist providing comfort.

### **Efficient Workflow**

Eliminates need for OR personnel to hold tubing.


### Low Profile Front Nozzle

*Improves visibility of the surgical site.* 

**Product Ordering Information** 

Description	Catalog Number	Quantity
PenAdapt Ideal for wall suction applications. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	PA1025	25/Case
PenAdapt (Cut-Back) Ideal for wall suction applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	PA1025CB	25/Case
PenAdapt Ideal for surgical smoke evacuator applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010	10/Case
PenAdapt Ideal for surgical smoke evacuator applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010CB	10/Case
PenAdapt Ideal for attachment to reusable electrosurgical pencil applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010B	10/Case
PenAdapt Ideal for attachment to reusable electrosurgical pencil applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010BCB	10/Case

<sup>\*</sup>ESU pencils sold separately.



PenAdape shown with EEU (reutsable) penci



PA2010B, PA2010BCB

