



PadPro®

Adult and Pediatric Multifunction
Electrodes for Defibrillation/Pacing/
Cardioversion/Monitoring




- Constructed to adhere and conform to patient's body, providing superior current distribution
- Adult electrodes offer a maximum of 50 shocks, 8 hours of continuous pacing or 24 hours of patient monitoring
- Pre-connect packaging options available
- 36-month shelf life
- PadPro® adapters, allowing multiple brands of defibrillators to use the same electrodes



ADULT
ECG ELECTRODES



PEDIATRIC
ECG ELECTRODES



PULSE OXIMETRY
SENSORS



To learn more about this and
other innovative products,
call **800-448-6506**
visit **Conmed.com** or email
info@mail.Conmed.com

Ordering Information

Adult and Pediatric Multifunction Electrodes for Defibrillation/Pacing/Cardioversion/Monitoring



Size A



Size A/B



Size C



Size D



Pre-Connect Packaging

Connector Style and Product Label Color

Product Description	Pouch Color	Background Label Color	Size	Universal/Anderson ¹	Zoll ²	Physio-Control [®] Quik-Combo™	Philips Barrel-Style	Philips Plug-Style	R2 ^{® 2}
Radiotransparent, Adult	Silver	White	A	2001	2001Z	2001M	2001P	2001H	2001R
Radiotransparent, Adult, w/carbon leadwires	Silver	White	A	2001-C	2001Z-C	2001M-C	2001P-C	2001H-C	N/A
Radiotransparent, Adult, w/Pre-Connect Packaging	Silver	White	A	2001-PC	2001Z-PC	2001M-PC	2001P-PC	2001H-PC	N/A
Radiotransparent, Adult, w/large-posterior pad	Silver	White	A/B	2001-EPS	N/A	N/A	N/A	N/A	N/A
Radiotransparent, Adult w/Sterile Packaging	Silver	White	A	2502	2502Z	2502M	2502P	N/A	N/A
Radiotranslucent, Adult	White	White	A	2516	2516Z	2516M	2516P	2516H	2516R
Radiotranslucent, Adult, w/Pre-Connect Packaging	White	White	A	2516-PC	2516Z-PC	2516M-PC	2516P-PC	2516H-PC	N/A
Radiotranslucent, Pediatric	White	Pink	D	2603	2603Z	2603M	2603P	2603H	2603R
Radiotranslucent, Mini Pediatric	White	Blue	C	2602	2602Z	2602M	2602P	2602H	N/A

Sizes:

A = 6.0" (15.2cm) x 4.25" (10.8cm)

B = 9.0" (22.9cm) x 6.0" (15.2cm)

C = 2.5" (6.4cm) x 2.125" (5.4cm)

D = 4.25" (10.8cm) x 2.875" (7.3cm)

¹ Electrodes with the Universal/Anderson connector can be utilized in environments with multiple brands of defibrillators (via PadPro adapters or PadPro cables).

² Electrodes with the R2 connector can be utilized in environments with multiple brands of defibrillators (via R2 adapters or R2 cables).

Product and company names listed are trademarks or trade names of their respective companies.