




Resting ECG Electrodes

With a variety of shapes, sizes, types, gels, and materials, ConMed has a resting ECG electrode for every monitoring application. We're one company with one goal – to satisfy the customer in the most cost effective method possible.



ADULT
ECG ELECTRODES



PEDIATRIC
ECG ELECTRODES



SPECIALTY
ECG ELECTRODES



To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.

RESTING ECG
ELECTRODES



ALLIGATOR CLIP
ADAPTER

Fastrace[®] 4 1915-100



Gel 2000 formulation sticks better and stays moist longer – even when the pouch has been opened.

- Silver/silver chloride design provides stable baselines and excellent electrical characteristics
- Conductive adhesive gel maintains contact with skin, even on hairy areas
- PVC free
- 24 month shelf life

V-Trace[®] 1935-100



Gel 2000 formulation sticks better and stays moist longer – even when the pouch has been opened.

- Silver/silver chloride design provides stable baselines and excellent electrical characteristics
- Conductive adhesive gel maintains contact with skin, even on hairy areas
- PVC free
- 24 month shelf life

HeartBeat[®] 1925-100



Gel 2000MP formulation provides improved staying power, but gentle release.

- Smaller design is ideal for adults and pediatrics
- Silver/silver chloride design provides stable baselines and excellent electrical characteristics
- Pliable material is more comfortable and prevents lift off
- PVC free
- 24 month shelf life

Snaptrace®

1920-050



The small circular design is ideal for adults and pediatrics, and offers easy placement for V-leads.

- Silver/silver chloride design provides stable baselines and excellent electrical characteristics
- Wet gel formulation for quick signal acquisition
- Radiolucent and MR safe
- PVC free
- 24 month shelf life

Snaptrace® II

1921-100



The smaller design is ideal for adults and pediatrics, and offers easy placement for V-leads.

- Silver/silver chloride design provides stable baselines and excellent electrical characteristics
- Solid adhesive hydrogel provides staying power but gentle release
- Radiolucent and MR safe
- PVC free
- 24 month shelf life

RESTING ECG
ELECTRODES



UNIVERSAL ADAPTER

Ordering Information

Ref. No.	Description	Electrodes/Pouch	Pouches/Case	Electrodes/Case
Tab Style Electrodes				
1915-100	Fastrace® 4	100	10	1,000
1935-100	V-Trace®	100	40	4,000
1925-100	Heartbeat	100	10	1,000
Snap Style Electrodes				
1920-050	Snaptrace®	50	20	1,000
1921-100	Snaptrace® II	100	10	1,000

Ref. No.	Fits Banana .156" (4.0mm)	Fits Pin .117" (3.0mm)	Fits Pin .125" (3.2mm)	Qty/Pkg	Description	
Tab Style Electrode Adapters						
ACB	●			10	Alligator clip adapter 6" (15.2cm) long	
ACP			●	10		
HPB-B	●			10	Bidirectional clip adapter 6" (15.2cm) long	
HPP-B			●	10		
HPH-Blue	●			10	One piece adapter	
HPH-Red		●	●	10		
Snap Style Electrode Adapters						
RLB-10	●			10	Grabber adapters 6" (15.2cm) long	
RLP-10			●	10		
SNB	●			10		
SNF		●		10	Snap adapters 6" (15.2cm) long	
SNP			●	10		
Universal Adapters for Tab or Snap Style Electrodes						
UMCL	●	●	●	10	Universal adapter	

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

ConMed.com
info@mail.conmed.com