

HELPING YOU PREVENT LOWER RESPIRATORY TRACT INFECTIONS



SOLUTIONS FOR RESPIRATORY HEALTH

Avanos is committed to help you reduce lower respiratory tract infections by providing effective clinical solutions to help with the removal and prevention of aspiration and colonisation in the oropharynx.

Avanos's Respiratory Health solutions bundle includes the BALLARD[®] Closed Suction Systems, MICROCUFF[®] Endotracheal Tubes, BALLARD[®] Oral Care Systems and BALLARD[®] Mini-BAL Sampling Catheter.




TABLE OF CONTENTS

CLOSED SUCTION SYSTEMS - ADULTS	5
BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS	5
BALLARD* Turbo-Cleaning Closed Suction System - Double Swivel Elbows	5
BALLARD* Turbo-Cleaning Closed Suction System - WET PAK* Double Swivel Elbows	7
BALLARD* Turbo-Cleaning Closed Suction System - T-Pieces	8
BALLARD* Turbo-Cleaning Closed Suction System - WET PAK* T-Pieces	9
BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS	10
BALLARD* Multi-Access Port Closed Suction System - Double Swivel Elbows	10
BALLARD* Multi-Access Port Closed Suction System - WET PAK* Double Swivel Elbows	12
BALLARD* Multi-Access Port Closed Suction Replacement Catheter	14
BALLARD* CLOSED SUCTION SYSTEMS	15
BALLARD* Closed Suction System - T-Pieces	15
BALLARD* Closed Suction System - Double Swivel Elbows	16
BALLARD* Closed Suction System - WET PAK* Double Swivel Elbows	17
BALLARD* CLOSED SUCTION SYSTEMS - ACCESSORIES	18
BALLARD* Multi-Access Port Closed Suction System - Accessories	18
BALLARD* Closed Suction System - Accessories	18
BALLARD* Closed Ventilator Circuit Drain	19
ENDOTRACHEAL TUBES - ADULTS	20
MICROCUFF* Subglottic Suctioning Endotracheal Tubes	20
MICROCUFF* Endotracheal Tubes	20
ORAL CARE SOLUTIONS - ADULTS	21
BALLARD* ORAL CARE KITS	21
BALLARD* Oral Care Advanced Kit	21
BALLARD* Oral Care Moderate Kit	21
BALLARD* Mini Oral Care Kit	22
BALLARD* ORAL CARE SUCTION KITS - COMPONENTS & SUCTION KITS	23
BALLARD* Oral Care Suction Components	23
BALLARD* Oral Care Swabs	24
BALLARD* Mini Oral Care Swabs	24
DIAGNOSTIC CATHETERS - ADULTS	25
BALLARD* MINI-BAL SAMPLING CATHETER	25

CLOSED SUCTION SYSTEMS - NEONATES/PAEDIATRICS	26
BALLARD* Closed Suction System for Neonates/Paediatrics - Elbows	26
BALLARD* Closed Suction System for Neonates/Paediatrics - "Y"	27
BALLARD* Multi-Access Catheter for Neonates/Paediatrics	28
BALLARD* Closed Suction System for Neonates/Paediatrics - Accessories	28
ENDOTRACHEAL TUBES - PAEDIATRICS	29
MICROCUFF* Endotracheal Tube for Paediatric, Oral/Nasal Magill	29
MICROCUFF* Endotracheal Tube for Paediatric, Oral Curved	30
ORAL CARE SOLUTIONS - NEONATES/PAEDIATRICS	31
BALLARD* Oral Care Suction Catheters for Neonates/Paediatrics	31
INDEX BY CODE	32


CLOSED SUCTION SYSTEMS - ADULTS

BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS

Indication: Aspiration catheter for removal of tracheal secretions of intubated / tracheotomized adult patients respecting the closed ventilation circuit.

BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEM - DOUBLE SWIVEL ELBOWS

- BALLARD* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 5 Years Shelf Life
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Available with dual lumen catheter
- Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex
- Does not contain DEHP
- Connector 15mm Female to 22mm Female



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
DOUBLE SWIVEL ELBOW	2271018-5	10 Fr/3.3mm	54cm	black	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	227101-5	10 Fr/3.3mm	54cm	black	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	22716-5	12 Fr/4mm	54cm	white	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	2271603-5	12 Fr/4mm	30.5cm	white	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	2271218-5	12 Fr/4mm	54cm	white	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	22712183	12Fr/4mm	30.5cm	white	20/case
Tracheostomy MDI Port					
DOUBLE SWIVEL ELBOW	227-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal					


All codes in **bold** indicate a sterile product

DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
DOUBLE SWIVEL ELBOW	22703-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	2271418-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	22714183-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy MDI Port					
DOUBLE SWIVEL ELBOW	227161-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	2271613-5	16 Fr/5.3mm	30.5cm	orange	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	2271618-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	22716183-5	16 Fr/5.3mm	30.5cm	orange	20/case
Tracheostomy MDI Port					

All codes in **bold** indicate a sterile product

BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEM - WET PAK* DOUBLE SWIVEL ELBOWS

- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Does not contain Natural Rubber Latex nor DEHP
- Validated for a maximum Usage Time of 72 Hours




DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
WET PAK* DOUBLE SWIVEL ELBOW	227166-5	12 Fr/4mm	54cm	white	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	22716036-5	12 Fr/4mm	30.5cm	white	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	22712186-5	12 Fr/4mm	54cm	white	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	227121863-5	12 Fr/4mm	30.5cm	white	20/case
Tracheostomy MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	2276-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	227036-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	22714186-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	227141863-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	2271616-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	22716186-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	227161863-5	16 Fr/5.3mm	30.5cm	orange	20/case
Tracheostomy MDI Port					

All codes in **bold** indicate a sterile product

BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEM - T-PIECES

- Gamma Sterilised: 5 Years Shelf Life
- Does not contain Natural Rubber Latex nor DEHP
- Validated for a maximum Usage Time of 72 Hours
- 22mm x 22mm Flex Tube Connector
- 15mm x 22mm Swivel Adapter (except 227105-5°)




DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
T-PIECE Endotracheal	227105-5°	10 Fr/3.3mm	54cm	black	20/case
T-PIECE Tracheostomy	2271053-5	10 Fr/3.3mm	30.5cm	black	20/case
T-PIECE Tracheostomy	227135-5	12 Fr/4mm	30.5cm	white	20/case
T-PIECE Endotracheal	22715-5	12Fr/4mm	54cm	white	20/case
T-PIECE Tracheostomy	2270135	14Fr/4.6mm	30.5cm	green	20/case

All codes in **bold** indicate a sterile product

BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEM - WET PAK* T-PIECES

- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
WET PAK* T-PIECE	2271516-5	12 Fr/4mm	54cm	white	20/case
Endotracheal					
WET PAK* T-PIECE	2270156-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal					
WET PAK* T-PIECE	22701356-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
WET PAK* T-PIECE	22012356-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					


All codes in **bold** indicate a sterile product

BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS

Indication: Aspiration catheter for intubated or tracheotomized adult patients who are on mechanical ventilation allowing doctors to perform alternate respiratory treatments without breaching the circuit.

BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM - DOUBLE SWIVEL ELBOWS

- BALLARD* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 5 Years Shelf Life
- Rotational Manifold with Multiple Ports (Suctioning Port & Alternate Therapy for Bronchoscopy or BALLARD* Mini-BAL sampling catheter)
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
DOUBLE SWIVEL ELBOW	24010-5	10 Fr/3.3mm	58.5cm	black	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	240121-5	12 Fr/4mm	58.5cm	white	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	240128-5	12 Fr/4mm	58.5cm	white	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	240123-5	12 Fr/4mm	35cm	white	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	24001238-5	12 Fr/4mm	35cm	white	20/case
Tracheostomy MDI Port					
DOUBLE SWIVEL ELBOW	24014-5	14 Fr/4.6mm	58.5cm	green	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	240148-5	14 Fr/4.6mm	58.5cm	green	20/case
Endotracheal MDI Port					


All codes in **bold** indicate a sterile product

DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
DOUBLE SWIVEL ELBOW	240143-5	14 Fr/4mm	35cm	green	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	2401438-5	14 Fr/4.6mm	35cm	green	20/case
Tracheostomy MDI Port					
DOUBLE SWIVEL ELBOW	240168-5	16 Fr/5.3mm	58.5cm	orange	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	2401638-5	16 Fr/5.3mm	35cm	orange	20/case
Tracheostomy MDI Port					

All codes in **bold** indicate a sterile product

BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM - WET PAK* DOUBLE SWIVEL ELBOWS

- BALLARD* Technology
- Built-in Cleaning Chamber
- Rotational Manifold with Multiple Ports (Suctioning Port & Alternate Therapy for Bronchoscopy or BALLARD* Mini-BAL sampling catheter)
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
WET PAK* DOUBLE SWIVEL ELBOW	240126-5	12 Fr/4mm	58.5cm	white	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	2401236-5	12 Fr/4mm	35cm	white	20/case
Tracheostomy					
WET PAK DOUBLE SWIVEL ELBOW	2401286-5	12 Fr/4mm	58.5cm	white	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	240146-5	14 Fr/4.6mm	58.5cm	white	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	2401436-5	14 Fr/4.6mm	35cm	green	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	2401486-5	14 Fr/4.6mm	58.5cm	green	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	24014386-5	14 Fr/4.6mm	35cm	green	20/case
Tracheostomy MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	240166-5	16 Fr/5.3mm	58.5cm	orange	20/case
Endotracheal					

All codes in **bold** indicate a sterile product

DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
WET PAK* DOUBLE SWIVEL ELBOW	2401636-5	16 Fr/5.3mm	35cm	orange	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	2401686-5	16 Fr/5.3mm	58.5cm	orange	20/case
Endotracheal MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	24016386-5	16 Fr/5.3mm	35cm	orange	20/case
Tracheostomy MDI Port					

All codes in **bold** indicate a sterile product

BALLARD* MULTI-ACCESS PORT CLOSED SUCTION REPLACEMENT CATHETER

- BALLARD* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 3 Years Shelf Life
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
ENDOTRACHEAL Replacement Catheter	240125-5	12 Fr/4mm	58.5cm	white	20/case
TRACHEOSTOMY Replacement Catheter	2401253-5	12 Fr/4mm	35cm	white	20/case
ENDOTRACHEAL Replacement Catheter	240145-5	14 Fr/4.6mm	58.5cm	green	20/case
TRACHEOSTOMY Replacement Catheter	2401453-5	14 Fr/4.6mm	35cm	green	20/case
ENDOTRACHEAL Replacement Catheter	240165-5	16 Fr/5.3mm	58.5cm	orange	20/case


All codes in **bold** indicate a sterile product

BALLARD* CLOSED SUCTION SYSTEMS

Indication: Aspiration catheter for removal of tracheal secretions of intubated / tracheotomized adult patients respecting the closed ventilation circuit.

BALLARD* CLOSED SUCTION SYSTEM - T-PIECES

- BALLARD* Technology
- Tactile Sleeve for Better Feel
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Validated for a maximum Usage Time of 24 Hours
- 22mm x 22mm Flex Tube Connector
- Available for Endotracheal and Tracheostomy tubes
- Available with Directional Tip
- Available with Dual Lumen Catheter
- Does not contain Natural Rubber Latex nor DEHP




DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
CLEAR T-PIECE Endotracheal	2205-5	14 Fr/4.6mm	54cm	green	20/case
CLEAR T-PIECE Endotracheal Directional Tip	22059-5	14 Fr/4.6mm	56cm	green	20/case
CLEAR T-PIECE Tracheostomy 15mm x 22mm Swivel Adapter	220135-5	14 Fr/4.6mm	30.5cm	green	20/case

All codes in **bold** indicate a sterile product

BALLARD* CLOSED SUCTION SYSTEM - DOUBLE SWIVEL ELBOWS

- BALLARD* Technology
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Gamma Sterilised: 5 Years Shelf Life
- Available for Endotracheal and Tracheostomy tubes
- Available with Dual Lumen Catheter
- Available with Directional Tip
- Connector 15mm Female to 22mm Female
- Validated for a maximum Usage Time of 24 Hours
- Does not contain Natural Rubber Latex nor DEHP




DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
DOUBLE SWIVEL ELBOW	2160-5	12 Fr/4mm	54cm	white	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	21609-5	12 Fr/4mm	56cm	white	20/case
Endotracheal Directional Tip					
DOUBLE SWIVEL ELBOW	21603-5	12 Fr/4mm	30.5cm	white	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	2210-5	14 Fr/4.6mm	54cm	white	20/case
Endotracheal					
DOUBLE SWIVEL ELBOW	22103-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
DOUBLE SWIVEL ELBOW	22107-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal Dual Lumen					
DOUBLE SWIVEL ELBOW	22108-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal MDI Port					
DOUBLE SWIVEL ELBOW	22109-5	14 Fr/4.6mm	56cm	green	20/case
Endotracheal Directional Tip					
DOUBLE SWIVEL ELBOW	2260-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal					

All codes in **bold** indicate a sterile product

BALLARD* CLOSED SUCTION SYSTEM - WET PAK* DOUBLE SWIVEL ELBOWS

- BALLARD* Technology
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Tactile Sleeve for Better Feel
- Contains in addition 12 Saline Vials
- Validated for a maximum Usage Time of 24 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
WET PAK* DOUBLE SWIVEL ELBOW	21606-5	12 Fr/4mm	54cm	white	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	216036-5	12 Fr/4mm	30.5cm	white	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	22106-5	14 Fr/4.6mm	54cm	green	20/case
Endotracheal					
WET PAK* DOUBLE SWIVEL ELBOW	221036-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
WET PAK* DOUBLE SWIVEL ELBOW	221069-5	14 Fr/4.6mm	56cm	green	20/case
Endotracheal					
Directional Tip					
WET PAK* DOUBLE SWIVEL ELBOW	2210386-5	14 Fr/4.6mm	30.5cm	green	20/case
Tracheostomy					
MDI Port					
WET PAK* DOUBLE SWIVEL ELBOW	22606-5	16 Fr/5.3mm	54cm	orange	20/case
Endotracheal					

All codes in **bold** indicate a sterile product

BALLARD* CLOSED SUCTION SYSTEMS - ACCESSORIES

BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM - ACCESSORIES


- Gamma Sterilised




DESCRIPTION	CODE	PACKAGING
BALLARD* SEAL CASSETTE	24013-5	20/case

BALLARD* CLOSED SUCTION SYSTEM - ACCESSORIES

- Gamma Sterilised



118




119

DESCRIPTION	CODE	PACKAGING
FLEX TUBE, 22MM X 22MM	113	50/case
Y-ADAPTER FOR SUCTION CANISTER	114	50/case
UNIT DOSE SALINE, 0.9% SODIUM CHLORIDE SOLUTION, PINK VIAL, 15ML, NONPYROGENIC, ASEPTICALLY FILLED	118	144/case
CONNECTOR 15MM FEMALE TO 22MM FEMALE	1115	50/case
SWIVEL ADAPTER, 15MM X 22MM, T-PIECE	111	50/case
FIXED ADAPTER, 15MM X 22MM, T-PIECE	112	50/case
METERED DOSE INHALER ADAPTER	119	20/case
BLUE ADAPTER FOR BALLARD* CLOSED SUCTION SYSTEM	50040	20/case

All codes in **bold** indicate a sterile product

BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN

- Allows the Clearing of Ventilator Condensation Without Exposing the Caregiver or the Patient
- Vacuum Activated
- Maintains a Closed System
- Developed With the Safety of the Caregiver and Patient in Mind



DESCRIPTION	CODE	PACKAGING
BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN ONLY	131	50/case
BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN KIT:	132	25/case
Circuit Drain 2,74 m Connector Tubing Y-Adapter		
BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN PROBE	133	50/case
2,74 m Connector Tubing Y-Adapter		
BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN WITH 90 DEGREE CONNECTION	134	25/case
2,74 m Connector Tubing Collapsible Connector Y-Adapter		

All codes in **bold** indicate a sterile product


ENDOTRACHEAL TUBES - ADULTS

Indication: Airway management by nasal or oral intubation of the trachea in adult patients

- Oral/Nasal Magill
- Helps to prevent micro-aspiration
- Sealed airway at very low cuff pressure
- EO Sterilised: 5 Years Shelf Life
- Radio-opaque
- Anatomically-based Intubation Depth Mark
- Murphy Eye
- Does not contain Natural Rubber Latex

MICROCUFF* SUBGLOTTIC SUCTIONING ENDOTRACHEAL TUBES

- Available in sizes 7.0 to 9.0mm (inner diameter)
- Cylindrical polyurethane ultra-thin cuff
- DEHP free



TUBE SIZE (INNER DIAMETER)	CODE	PACKAGING
7.0mm	13220	1 case - 10 eaches/dispenser
7.5mm	13221	1 case - 10 eaches/dispenser
8.0mm	13222	1 case - 10 eaches/dispenser
8.5mm	13223	1 case - 10 eaches/dispenser
9.0mm	13224	1 case - 10 eaches/dispenser

MICROCUFF* ENDOTRACHEAL TUBES

- Available in sizes 5.0 to 10.0mm (inner diameter)
- Cylindrical Cuff



TUBE SIZE (INNER DIAMETER)	CODE	PACKAGING
5.0mm	35210	1 case - 10 eaches/dispenser
5.5mm	35211	1 case - 10 eaches/dispenser
6.0mm	35212	1 case - 10 eaches/dispenser
6.5mm	35213	1 case - 10 eaches/dispenser
7.0mm	35214	1 case - 10 eaches/dispenser
7.5mm	35215	1 case - 10 eaches/dispenser
8.0mm	35216	1 case - 10 eaches/dispenser
8.5mm	35217	1 case - 10 eaches/dispenser
9.0mm	35218	1 case - 10 eaches/dispenser
10.0mm	35220	1 case - 10 eaches/dispenser

All codes in **bold** indicate a sterile product

ORAL CARE SOLUTIONS - ADULTS

BALLARD* ORAL CARE KITS

Indication: Suction tools for removal of oral secretions of patients, that support the compliance to a comprehensive oral care protocol.



BALLARD* ORAL CARE ADVANCED KIT

- 2 Years Shelf Life

Kit includes:

- 1 Prep pack (incl. 1 covered Yankauer, 1 Y-adapter & 1 suction handle)
- 4 Suction Swab Packs (incl. 1 suction swab, 1 applicator swab & 1 empty solution cup)
- 2 Toothbrush Packs (incl. 1 suction toothbrush, 1 oral swab & 1 empty solution cup)

DESCRIPTION	CODE	PACKAGING
BALLARD* ORAL CARE ADVANCED KIT	25928	24/case

BALLARD* ORAL CARE MODERATE KIT

- 2 Years Shelf Life

Kit includes:

- 1 Prep Pack (incl. 1 covered Yankauer, 1 Y-adapter & 1 suction handle)
- 2 Suction Swab Packs (incl. 1 suction swab, 1 applicator swab & 1 empty solution cup)
- 1 Toothbrush Pack (incl. 1 suction toothbrush, 1 oral swab & 1 empty solution cup)

DESCRIPTION	CODE	PACKAGING
BALLARD* ORAL CARE MODERATE KIT	25930	32/case

All codes in **bold** indicate a sterile product

BALLARD* MINI ORAL CARE KIT

- 2 Years Shelf Life
- For use with a 4-hour oral care protocol



Kit includes:

- 1 Prep Pack (reusable for 24 Hours - incl. self-cleaning covered Yankauer)
- 2 Toothbrush Packs (incl. 1 suction toothbrush, 1 mini applicator swab, 1 empty solution cup & 1 mini suction catheter)
- 4 Suction Swab Packs (incl. 1 mini suction swab, 1 mini applicator swab & 1 empty solution cup)

DESCRIPTION	CODE	PACKAGING
BALLARD* MINI ORAL CARE KIT	970008-5	24/case


All codes in **bold** indicate a sterile product

BALLARD* ORAL CARE SUCTION KITS - COMPONENTS & SUCTION KITS

Indication: Suction tools for removal of oral secretions of patients, that support the compliance to a comprehensive oral care protocol.

BALLARD* ORAL CARE SUCTION COMPONENTS

- 5 Years Shelf Life (except 99785[◊], 99786[†], 99789[‡]: 2 Years Shelf Life)




DESCRIPTION	CODE	PACKAGING
BALLARD* READY CARE* ORAL CARE SUCTION PROBES	123	50/case
BALLARD* ORAL CARE SUCTION CATHETER FOR ADULTS	12601	25/case
BALLARD* ORAL CARE SUCTION TOOTHBRUSH	12602	25/case
BALLARD* ORAL CARE SUCTION SWAB, ANGLED-TIP	12603	25/case
BALLARD* SUCTION VALVE AND Y-ADAPTER	12606	25/case
BALLARD* ORAL CARE SUCTION SWAB, WITH DENTIFRICE (MAY BE USED WITHOUT BALLARD* SUCTION VALVE)	12607	25/case
BALLARD* SELF-CLEANING COVERED YANKAUER	99785 [◊]	25/case
BALLARD* SUCTION HANDLE	99786 [†]	25/case
BALLARD* ORAL CARE PREP PACK:	99789 [‡]	40/case
1 self-cleaning, covered Yankauer, 1 suction handle, Y-adapter		

All codes in **bold** indicate a sterile product

BALLARD* ORAL CARE SWABS

- Shelf Life for Swabs with Dentifrice: 2 years
- Shelf Life for Swabs without Dentifrice: 5 years



DESCRIPTION	CODE	ITEMS/POUCH	POUCHES/CASE
SINGLE WRAPPED			
BALLARD* ORAL CARE SWAB	12240	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	12245	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	12246	250	2
BALLARD* ORAL CARE SWAB, PLASTIC HANDLE	12254	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE, PLASTIC HANDLE	12255	250	4
BULK PACKAGED			
BALLARD* ORAL CARE SWAB	12241	500	2
BALLARD* ORAL CARE SWAB	12243	20	50
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	12249	20	50
BALLARD* ORAL CARE SWAB, PLASTIC HANDLE	12251	20	50
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE, PLASTIC HANDLE	12257	20	50

BALLARD* MINI ORAL CARE SWABS

- Shelf Life for Swabs: 5 years
- Individually Wrapped with 12cm Plastic Handle

DESCRIPTION	CODE	ITEMS/POUCH	POUCHES/CASE
BALLARD* MINI ORAL CARE SWABS	122548	250	4


All codes in **bold** indicate a sterile product

DIAGNOSTIC CATHETERS - ADULTS

BALLARD* MINI-BAL SAMPLING CATHETER

Indication: This product is intended for early, accurate diagnosis of VAP.

- Simple, effective non-bronchoscopic bronchoalveolar lavage at the bedside
- EO Sterilised: 3 Years Shelf Life
- Directional Tip for R/L Insertion
- Includes an access port elbow, adapter for seal cassette, seal cassette, stepped connector, oxygen line
- Atraumatic "Mushroom" Tip
- Does not contain Natural Rubber Latex



DESCRIPTION	CODE	PACKAGING
BALLARD* MINI-BAL SAMPLING CATHETER - 16 FR OUTER/12 FR INNER TELESCOPING CATHETER	142	5/case
BALLARD* MINI-BAL SAMPLING CATHETER - 13 FR OUTER/8 FR INNER TELESCOPING CATHETER	143	5/case


All codes in **bold** indicate a sterile product

CLOSED SUCTION SYSTEMS - NEONATES/PAEDIATRICS

Indication: Aspiration catheter for removal of tracheal secretions of intubated neonatal/paediatric patients respecting the closed ventilation circuit.

BALLARD* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - ELBOWS

- BALLARD* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Contains Elbow Adapters
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
ELBOW ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm)	206-5	6 Fr/2mm	30.5cm	light green	20/case
ELBOW ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm, 4mm)	207-5	7 Fr/2.3mm	35.5cm	ivory	20/case
ELBOW ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	208-5	8 Fr/2.6mm	35.5cm	light blue	20/case
ELBOW ENDOTRACHEAL TUBE ADAPTERS (4.5mm, 5mm, 5.5mm)	210-5	10 Fr/3.3mm	40.5cm	black	20/case

All codes in **bold** indicate a sterile product

BALLARD* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - "Y"

- BALLARD* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Contains Elbow and Y-adapters
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
"Y" ENDOTRACHEAL TUBE ADAPTERS (2mm, 2.5mm)	195-5	5 Fr/1.6mm	30.5cm	grey	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm)	196-5	6 Fr/2mm	30.5cm	light green	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	197-5	7 Fr/2.3mm	30.5cm	ivory	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	198-5	8 Fr/2.6mm	30.5cm	light blue	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (4.5mm, 5mm, 5.5mm)	1910-5	10 Fr/3.3mm	40.5cm	black	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (5.5mm, 6mm, 6.5mm)	1912-5	12 Fr/4mm	40.5cm	white	20/case

All codes in **bold** indicate a sterile product

BALLARD* MULTI-ACCESS CATHETER FOR NEONATES/PAEDIATRICS


Indication: A multi-access catheter designed to deliver surfactants to the distal airway.

- BALLARD* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Contains Closure Cap and Y-adapter
- Does not contain Natural Rubber Latex
- The Catheter is DEHP free



DESCRIPTION	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
MAC CATHETER ADAPTER (2.5mm, 3mm, 3.5mm)	1900	5 Fr/1.6mm	30.5cm	white	5/case
MAC CATHETER ADAPTER (2mm)	1920	5 Fr/1.6mm	30.5cm	white	5/case
MAC CATHETER ADAPTER (4mm)	1940	5 Fr/1.6mm	30.5cm	white	5/case
TRANSPORT KIT ("Y" Plug and Cap)	19000				50/case

BALLARD* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - ACCESSORIES



DESCRIPTION	CODE	PACKAGING
ENDOTRACHEAL TUBE ADAPTER 2MM WITH LAVAGE PORT	2-020	20/case
ENDOTRACHEAL TUBE ADAPTER 2.5MM WITH LAVAGE PORT	2-025	20/case
ENDOTRACHEAL TUBE ADAPTER 3MM WITH LAVAGE PORT	2-030	20/case
ENDOTRACHEAL TUBE ADAPTER 3.5MM WITH LAVAGE PORT	2-035	20/case
ENDOTRACHEAL TUBE ADAPTER 4MM WITH LAVAGE PORT	2-040	20/case
ENDOTRACHEAL TUBE ADAPTER 4.5MM WITH LAVAGE PORT	2-045	20/case
MANUAL PERCUSSOR FOR NEONATAL PULMONARY USE	11000	12/case


All codes in **bold** indicate a sterile product

ENDOTRACHEAL TUBES - PAEDIATRICS

Indication: Airway management by nasal or oral intubation of the trachea in paediatric patients.

- Available in Magill & RAE / preformed types
- Available in sizes 3.0 to 7.0mm (inner diameter)
- EO Sterilised: 5 Years Shelf Life
- Designed for ideal anatomic cuff placement
- Sealed airway at ultra low cuff pressure
- Does not contain Natural Rubber Latex


MICROCUFF* ENDOTRACHEAL TUBE FOR PAEDIATRIC, ORAL/NASAL MAGILL



TUBE SIZE (INNER DIAMETER)	CODE	AGE/WEIGHT YEARS/KG	PACKAGING
3.0mm	35111	Term $\geq$ 3kg to < 8 months	10/dispenser
3.5mm	35112	8 months < 2 years	10/dispenser
4.0mm	35113	2 to < 4 years	10/dispenser
4.5mm	35114	4 to < 6 years	10/dispenser
5.0mm	35115	6 to < 8 years	10/dispenser
5.5mm	35116	8 to < 10 years	10/dispenser
6.0mm	35117	10 to < 12 years	10/dispenser
6.5mm	35118	12 to < 14 years	10/dispenser
7.0mm	35119	14 to < 16 years	10/dispenser

All codes in **bold** indicate a sterile product

MICROCUFF* ENDOTRACHEAL TUBE FOR PAEDIATRIC, ORAL CURVED



TUBE SIZE (INNER DIAMETER)	CODE	AGE/WEIGHT YEARS/KG	PACKAGING
3.0mm	35161	Term ≥ 3kg to < 8 months	10/dispenser
3.5mm	35162	8 months < 2 years	10/dispenser
4.0mm	35163	2 to < 4 years	10/dispenser
4.5mm	35164	4 to < 6 years	10/dispenser
5.0mm	35165	6 to < 8 years	10/dispenser
5.5mm	35166	8 to < 10 years	10/dispenser
6.0mm	35167	10 to < 12 years	10/dispenser
6.5mm	35168	12 to < 14 years	10/dispenser
7.0mm	35169	14 to < 16 years	10/dispenser

All codes in **bold** indicate a sterile product

ORAL CARE SOLUTIONS - NEONATES/PAEDIATRICS

BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS

- Does not contain Natural Rubber Latex
- Packaging: 5 pouches of 50 units

DESCRIPTION	CODE	PACKAGING
5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS (6 Fr - 2mm)	1221	50/case
5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS (8 Fr - 2.6mm)	1222	50/case
5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS (10 Fr - 3.3mm)	1223	50/case

All codes in **bold** indicate a sterile product

INDEX BY CODE

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
2-020	28	12240	24	35116	29	240145-5	14
2-025	28	12241	24	35117	29	240146-5	12
2-030	28	12243	24	35118	29	240148-5	10
2-035	28	12245	24	35119	29	240165-5	14
2-040	28	12246	24	35161	30	240166-5	12
2-045	28	12249	24	35162	30	240168-5	11
35cm	10	12251	24	35163	30	970008-5	22
111	18	12254	24	35164	30	2210386-5	17
112	18	12255	24	35165	30	2270135	8
113	18	12257	24	35166	30	2270156-5	9
114	18	12601	23	35167	30	2271018-5	5
118	18	12602	23	35168	30	2271053-5	8
119	18	12603	23	35169	30	2271218-5	5
123	23	12606	23	35210	20	2271418-5	6
131	19	12607	23	35211	20	2271516-5	9
132	19	13220	20	35212	20	2271603-5	5
133	19	13221	20	35213	20	2271613-5	6
134	19	13222	20	35214	20	2271616-5	7
142	25	13223	20	35215	20	2271618-5	6
143	25	13224	20	35216	20	2401236-5	12
195-5	27	19000	28	35217	20	2401253-5	14
196-5	27	21603-5	16	35218	20	2401286-5	12
197-5	27	21606-5	17	35220	20	2401436-5	12
198-5	27	21609-5	16	50040	18	2401438-5	11
206-5	26	22059-5	15	99785◇	23	2401453-5	14
207-5	26	22103-5	16	99786†	23	2401486-5	12
208-5	26	22106-5	17	99789‡	23	2401636-5	13
210-5	26	22107-5	16	122548	24	2401638-5	11
227-5	5	22108-5	16	216036-5	17	2401686-5	13
1115	18	22109-5	16	220135-5	15	22012356-5	9
1221	31	22606-5	17	221036-5	17	22701356-5	9
1222	31	22703-5	6	221069-5	17	22712183	5
1223	31	22715-5	8	227036-5	7	22712186-5	7
1900	28	22716-5	5	227101-5	5	22714183-5	6
1910-5	27	24010-5	10	227105-5◇	8	22714186-5	7
1912-5	27	24013-5	18	227135-5	8	22716036-5	7
1920	28	24014-5	10	227161-5	6	22716183-5	6
1940	28	25928	21	227166-5	7	22716186-5	7
2160-5	16	25930	21	240121-5	10	24001238-5	10
2205-5	15	35111	29	240123-5	10	24014386-5	12
2210-5	16	35112	29	240125-5	14	24016386-5	13
2260-5	16	35113	29	240126-5	12	227121863-5	7
2276-5	7	35114	29	240128-5	10	227141863-5	7
11000	28	35115	29	240143-5	11	227161863-5	7

AVANOS

For more information, please call Avanos
at +32 2 700 68 51 - fax +32 2 711 26 91
or send an email to customerservice.export@avanos.com

www.avanos.co.uk