





Compatible Sensors

Pulse Oximetry Sensors are designed for use with a variety of pulse oximeter equipment. Finally there are choices with quality, performance and cost-savings.




ECG
ELECTRODES



BLOOD PRESSURE
CUFFS



CABLES
AND LEADWIRES



To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.

Nellcor* Compatible SpO₂ Sensors by ConMed

ADULT



Adult Disposable Sensors 2311, 3311-F, 3311-V

CROSS REFERENCE TO NELLCOR* PART NUMBERS: D-25, MAX-A

- Available in Vinyl Tape and Foam
- 9 pin male connector
- For use on adults > 30kg
- 18" cable length

CHILD



Pediatric Disposable Sensors 2312, 3312-F, 3312-V

CROSS REFERENCE TO NELLCOR* PART NUMBERS: D-20, MAX-P

- Available in Vinyl Tape and Foam
- 9 pin male connector
- For use on pediatrics 10-50kg
- 18" cable length

ADULT



Adult Reusable Sensors 2010, 2210

CROSS REFERENCE TO NELLCOR* PART NUMBERS: DS100A, D-YS, D-YSE

- 9 pin male connector
- 36" cable length

CHILD



Infant Disposable Sensors 2313, 3313-C

CROSS REFERENCE TO NELLCOR* PART NUMBERS: I-20, MAX-I

- For use on infants 10 < 20kg
- 36" cable length
- 9 pin male connector

CHILD



Neonatal Disposable Sensors 2314, 3314-C

CROSS REFERENCE TO NELLCOR* PART NUMBERS: N-25, MAX-N

- For use on neonatals < 10kg
- 36" cable length
- 9 pin male connector

9-PIN MALE
CONNECTOR



Nellcor* compatible SpO₂ sensors
feature a 9-pin male connector

NELLCOR* OXIMAX ADAPTER CABLE 2500GM

Allows ConMed Dolphin Nellcor* Sensors to be used with Nellcor* Oximax equipped monitors

- FDA 510K cleared
- 12" cable length



ConMed Dolphin® Compatible SpO₂ Sensors

SENSORS



2331 | 2332 | 2333



2030 | 2031



2321 | 2323



2020 | 2220

Datex* Compatible SpO₂ Sensors by ConMed

Datex Comp Sens Adult – Dolphin	2331
Datex Comp Sens Pediatric – Dolphin	2332
Datex Comp Sens Infant – Dolphin	2333
Reusable Finger Clip – Datex	2030
Reusable Finger Sensor – Datex	2331

Ohmeda* Compatible SpO₂ Sensors by ConMed

Adult / Pediatric Disposable – Ohmeda	2321
Infant / Neonatal Disposable – Ohmeda	2323
Reusable Finger Clip – Ohmeda	2020
Y-Multisite Reusable – Ohmeda	2220



2060



2260

Novamatrix* Compatible SpO₂ Sensors by ConMed

Reusable Finger – Novamatrix	2060
Reusable Y Multisite – Nova	2260



2270

Nonin* Compatible SpO₂ Sensors by ConMed

Reusable Finger Clip – Nonin	2270
------------------------------	------

CABLES



2411 | 2421



2412 | 2422



2413 | 2423



2414



2416 | 2426 DM



2417 | 2427



2418 | 2428

6' ConMed SpO₂ Extension / Adapter Cables

6' Extension Cable (9 Pin)-Nellcor*	2411
6' Extension Cable (9 Pin)-Ohmeda*	2412
6' Extension Cable (9 Pin)-Datex*	2413
6' Extension Cable-CSI*	2414
6' Extension Cable-Novamatrix*	2416
6' Extension Cable (9 Pin to HP)-Nellcor*	2417
6' Extension Cable (9 Pin to SL)-Nellcor*	2418

10' ConMed SpO₂ Extension / Adapter Cables

10' Extension Cable (9 Pin)-Nellcor*	2421
10' Extension Cable (9 Pin)-Ohmeda*	2422
10' Extension Cable (9 Pin)-Datex*	2423
10' Extension Cable-Novamatrix*	2426 DM
10' Extension Cable (9 Pin to HP)-Nellcor*	2427
10' Extension Cable (9 Pin to SL)-Nellcor*	2428

OTHER

ConMed SpO₂ Accessories

Adult / Neonatal Tape Wraps	22XX
Infant / Neonatal Disposable – Ohmeda	22XX-1
Neonatal Tape W/Y-Wraps	22XX-2
Ear Clip W/Y-Sensor	22XX-E



525 French Road, Utica, NY 13502

COMPATIBLE
SENSORS

Nellcor* Cross-Reference Information

Nellcor* Catalog Number	ConMed Catalog Number	Adapter Cable
D-25	2311, 3311-F, 3311-V	
D-20	2312, 3312-F, 3312-V	
I-20	2313, 3313-C	
N-25	2314, 3314-C	
MAX-A	2311, 3311-F, 3311-V	Use with 2500GM adapter cable
MAX-P	2312, 3312-F, 3312-V	Use with 2500GM adapter cable
MAX-I	2313, 3313-C	Use with 2500GM adapter cable
MAX-N	2314, 3314-C	Use with 2500GM adapter cable
DS100A	2010	Use with 2500GM adapter cable with Oximax* equipped monitors
D-YS	2010	Use with 2500GM adapter cable with Oximax* equipped monitors
D-YSE	2210	Use with 2500GM adapter cable with Oximax* equipped monitors

Other Manufacturer Cross-Reference Information

Manufacturer	Adult Reusable	Reusable Multisite Y	Adult Disposable	Pediatric Disposable	Infant Disposable	Neonatal Disposable
for Ohmeda*	2020	2220	2321	2321	2323	2323
for Datex*	2030	N/A	2331	2332	2333	N/A
for Novametrics*	2060	2260	N/A	N/A	N/A	N/A
for Nonin*	2070	N/A	N/A	N/A	N/A	N/A

* CSI, Dolphin, Nellcor, Ohmeda, Datex, Novametrics, Nonin and Oximax are trademarks of Cordis Systems, Inc., ConMed Corporation, Covidien LP, GE Healthcare Finland, RIC Investments, LLC, Nonin Medical, Inc., and Covidien LP, respectively.

Manufactured for ConMed by Spacelabs Healthcare Issaquah, WA.

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

ConMed.com
info@mail.conmed.com