

MEDIUM WEIGHT

Isolation Gown with Knit Cuffs



Skin
Friendly



Hydrophobic



Multi
Layer

Features/Benefits:

Dual Layer 30g Polypropylene Spunbond-Spunbond Fabric:

Protects against fluids and particles

Tear resistant

Breathable which prevents irritating
hot temperature build-up

Anti-Static Coating: Prevent the gown from clinging to your clothes

Front Waist Ties: Allow for easy and fast securement

Knit Cuffs

Size: One Size Fits Most, XL

Color: Blue, Green, Pink, White, Yellow

Quantity/UOM: 50pc/cs

Specifications

Protection Level	Created with water repellent
Material	30g Polypropylene Spunbond-Spunbond Fabric
Design	Neck & Front Waist Ties, Full Back
Cuff Style	Knit Cuffs
Size/Dimension	One Size Fits Most, XL - Length: 45" Width: 54"
Color	Blue, Green, Pink, White, Yellow
Quantity/UoM	50 Pieces/case

Item Number	Description
6010BC	MEDIUM WT. ISO GOWN KNIT CUFFS, BLUE
6010BC-XL	MEDIUM WT. ISO GOWN KNIT CUFFS, XL, BLUE
6010GC	MEDIUM WT. ISO GOWN KNIT CUFFS, GREEN
6010PC	MEDIUM WT. ISO GOWN KNIT CUFFS, PINK
6010WC	MEDIUM WT. ISO GOWN KNIT CUFFS, WHITE
6010YC	MEDIUM WT. ISO GOWN KNIT CUFFS, YELLOW



Knit Cuffs

