



Milliken



PRODUCT LINE OVERVIEW





THE COFLEX DIFFERENCE

For the Patient

- Improves comfort with lightweight and breathable material
- Conforms to area of application
- Stays in place

For the Application

- Secures primary dressings and devices on patients
- Compresses to help control swelling
- Protects patients with latex-free options

For the Clinician

- Improves clinical ease-of-use with EasyTear technology
- Delivers consistent application
- Offers HCPCS codes

Milliken's CoFlex cohesive bandages are designed to simplify and streamline the healthcare clinicians job.

USED IN:	USED FOR:	
Blood Banks, Labs, Donor Rooms, Plasmapheresis, Apheresis	Compression bandage	
Neonatal Intensive Care, Special Nurseries	Affix monitors to infants (1" Smiley prints are popular)	
Orthopedic Hand Therapy	Support and compression bandage	
Radiology, MRI, IV Therapy	Hold IV's in place	
Burn Centers	Pressure bandage to promote circulation and healing	
Occupational and Physical Therapy	Support bandage, edema control, and finger immobilization	
Operating Room	Vascular post-operative wrap securing drapes, tubes and equipment. Knee and shoulder arthroscopy, total knee post-operative dressing	
Dermatology	Compression/fixing bandage	
Orthopedics	Leg ulcer wrap with Unna boot as an overwrap	
Pediatrics	Compression bandage (Prints and colors are popular)	
Podiatry	Compression/fixing bandage	
Dialysis Units	Holding IV's and tubes	
Hemodialysis	Compression bandage	
Chiropractors	Compression bandage, heel spurs	
Plastic Surgery	Wrap to control swelling	







COFLEX LF²

Milliken's CoFlex LF² is ideal for high-risk patient populations and surgical patients. CoFlex LF² is designed to provide comfort to patients with sensitive skin, eliminate the risk associated with latex sensitivity and protect the sterile field in surgical procedures.

- Soft and comfortable gentle on the skin
- Ideal for OR application
- Water resistant for arthroscopic procedures**
- Sticks directly to a wet primary dressing
- EasyTear technology*
- Controlled compression will not constrict or loosen over time
- Stays in place does not shift or unravel with use



CoFlex LF² is available in single color cases or a color pack for variety!



Color Pack (CP) Includes: Neon Pink, Blue, Purple, Light Blue, Neon Green, Red

Catalog #	Size	Rolls/Case	Catalog #	Size	Rolls/Case
9100XX-030	1" x 5 yds.	30	9300(S)-024 sterile/tan	3" x 5 yds.	24
9100(S)-024 sterile/tan	1" x 5 yds.	24	9400XX-018	4" x 5 yds.	18
9150XX-048	1.5" x 5 yds.	48	9400(S)-018 sterile/tan	4" x 5 yds.	18
9200XX-036	2" x 5 yds.	36	9400L-018 (tan)	4" x 6.5 yds.	18
9200(S)-024 sterile/tan	2" x 5 yds.	24	9600XX-012	6" x 5 yds.	12
9300XX-024	3" x 5 yds.	24	9600(S)-012 sterile/tan	6" x 5 yds.	12







COFLEX NL

Milliken's CoFlex NL is a worry-free choice designed for use by everyone, everywhere. CoFlex NL provides excellent self adherence, excellent controlled compression, and excellent versatility to meet the multiple needs of an entire healthcare facility.

- Largest latex-free product offering with more than
 90 options
- EasyTear technology*
- Controlled compression will not constrict or loosen over time
- Stays in place does not shift or unravel with use
- Breathable and lightweight
- Made in the USA





CoFlex NL is available in single color cases or color packs for variety!

Color Pack (CP) Includes: Neon Pink, Blue, Purple, Light Blue, Neon Green, Red
Rainbow Pack (RB) Includes: Teal, Blue, White, Purple, Red, Green
Camo Color Pack (CMCP) Includes: Camo Print on Lt. Blue, Pink, Orange, & Green

Kids Pack (KP) Includes: Hearts, Dinosaurs, Cars & Trucks, Bumble Bees



Smiley Pack (SC) Includes: Smiley on Lavender, Lt. Blue, Neon Pink, & Yellow



Catalog #	Size	Rolls/Case	Catalog #	Size	Rolls/Case
5100XX-030	1" x 5 yds.	30	5400XX-018	4" x 5 yds.	18
5150XX-048	1.5" x 5 yds.	48	5400(S)-018 sterile/tan	4" x 5 yds.	18
5200XX-036	2" x 5 yds.	36	5400L-018 (tan)	4" x 6.5 yds.	18
5300XX-024	3" x 5 yds.	24	5600XX-012	6" x 5 yds.	12
5300(S)-024 sterile/tan	3" x 5 yds.	24	5600(S)-012 sterile/tan	6" x 5 yds.	12







COFLEX & COFLEX MED

Milliken's CoFlex & CoFlex Med are an economical solution when price is paramount. CoFlex and CoFlex Med contain latex and provide exceptional self adherence while providing controlled compression.

- Controlled compression will not constrict or loosen over time
- Stays in place does not shift or unravel with use
- Versatile and can be used for a wide variety of applications
- Breathable and lightweight
- EasyTear technology*
- Water resistant for arthroscopic procedures**




CoFlex and CoFlex Med are available in single color cases or color packs for variety!



Color Pack (CP) Includes: Neon Pink, Blue, Purple, Light Blue, Neon Green, Red Rainbow Pack (RB) Includes: Teal, Blue, White, Purple, Red, Green

CoFlex Catalog #	CoFlex Med Catalog #	Size	Rolls/Case	CoFlex Catalog #	CoFlex Med Catalog #	Size	Rolls/Case
3100XX-030	7100XX-030	1" x 5 yds.	30	3400XX-018	7400XX-018	4" x 5 yds.	18
3150XX-048	7150XX-048	1.5" x 5 yds.	48	-	7400(S)-024 sterile/tan	4" x 5 yds.	18
3200XX-036	7200XX-036	2" x 5 yds.	36	3400L-018 (tan)	-	4" x 6.5 yds.	18
3300XX-024	7300XX-024	3" x 5 yds.	24	3600XX-012	7600XX-012	6" x 5 yds.	12
-	7300(S)-024 sterile/tan	3" x 5 yds.	24	-	7600(S)-012 sterile/tan	6" x 5 yds.	12

YOU CONTROL THE TENSION



- 1. No to low tension
- 2. Medium tensior
 - 3. High tension


APPLICATION AND REMOVAL TIPS







2. Upon completion of wrapping, grasp bandage between fingers and remove stretch - handtear* or cut the bandage to desired length. Relax the last 2" of bandage and then press down.



3. Seal the end of the bandage wrap with your finger, lightly meshing the top of the bandage into the bottom layer. This will provide long staying power.

Removal: CoFlex bandages may be removed by unwrapping or with bandage scissors.



HCPCS CODES (Healthcare Common Procedure Coding System Codes)

Product	Size	HCPCS Code
CoFlex Med CoFlex NL CoFlex LF ²	1" x 5 yds. 1.5" x 5 yds. 2" x 5 yds.	A6453
CoFlex Med CoFlex NL CoFlex LF ²	3" x 5 yds. 4" x 5 yds.	A6454
CoFlex Med CoFlex NL CoFlex LF ²	6" x 5 yds.	A6455

Centers for Medicare & Medicaid Services (CMS) has verified that the HCPCS codes listed above have been assigned to Andover Healthcare, Inc.

(Updated July 2015)

Andover Healthcare, Inc., a subsidiary of Milliken Healthcare Products, LLC 9 Fanaras Drive, Salisbury, MA 01952 USA · www.healthcare.milliken.com · 800.432.6686

CoFlex is a registered trademark of Andover Healthcare, Inc., a subsidiary of Milliken Healthcare Products, LLC.

Document #COHESIVEPOr1

 $^{^{\}star}$ Easy Tear applicable for CoFlex Med, CoFlex NL and CoFlex LF 2 . Scissors required for CoFlex.

^{**}Water-resistant applicable for CoFlex Med and CoFlex LF².