



Milliken



PRODUCT LINE OVERVIEW





THE POWERFLEX DIFFERENCE

For the Athlete

- Supports muscles and joints
- Prevents injury
- Provides post-injury stability

For the Application

- Eliminates need for prewrap and adhesive spray
- Increases speed of application
- Compatible as the foundation for taping systems

For the Athletic Trainer

- Streamlines ordering by reducing materials needed from 3 to 1
- Allows for modification to address each athlete's needs
- Delivers consistent application

Andover by Milliken's cohesive wraps and tapes are designed to empower athletes to perform at their best while simplifying and streamlining the athletic trainers job.

In addition to these products, we support athletic trainers with literature and videos, sharing expertise on taping methods and specific injury types. These resources include a Taping Methods Video Textbook, which outlines 35 step-by-step techniques in preventative and post-injury taping instruction by Ron O'Neil. Ron is a formerly certified athletic trainer for 37 years – 29 of which were spent with the National Football League (NFL).

Backed by an organization with more than 150 years of experience in textile design, chemical science and materials manufacturing, we strive to improve clinical outcomes and enhance athletic performance while also helping our customers to increase efficiency, reduce waste and optimize their budgets.

Try our easy-to-use calculator to discover what your annual cost savings could be by converting to Andover by Milliken taping systems!





WHICH TAPE OR WRAP SHOULD I CHOOSE?

BEST Overall

ULTIMATE TAPING SYSTEM

Combining the strength and support of **POWERFLEX + POWERPRO**



Are YOU using PowerFlex as a base?

- 23 lbs. tensile strength
- Sweat & water resistant
- Replaces the need for prewrap and adhesive spray

Combine it with PowerPro for ultimate support.

- Versatile adhesive tape
- Less than 1% stretch
- Cotton free won't absorb sweat or water

BEST Strength

POWERFLEX & POWERTAPE

(cohesive athletic tape)

BEST Budget

POWERFLEX & VICTORYTAPE

(cohesive ahtletic tape)

BEST Speed

POWERSPEED & POWERPRO

(30% stronger than PowerFlex - apply less)

BEST Reinforcement

POWERFAST

(Apply over any system for additional support without restricting movement)

Build-Your-Own Taping System:

 TOP LAYER (tape)

POWERPRO Synthetic Adhesive Sports Tape
POWERTAPE Cohesive Athletic Tape
VICTORYTAPE Cohesive Athletic Tape

SELF-ADHERING SPORTS WRAP



POWERFLEX and POWERFLEX LF

PowerFlex® empowers stronger athletes with its superior performance, wide variety of applications and its compatibility as the foundation of taping systems. This strong, fabric-based wrap tears easily by hand and sticks to itself for quick wrapping and longer staying power.

Use PowerFlex as a base to streamline ordering. Replacing prewrap, adhesive spray, and heel & lace pads, it reduces materials needed from 3 to 1!

- 23 lbs. tensile strength
- Sweat & water resistant
- Eliminates need for prewrap and adhesive spray
- Supports muscles & joints
- Provides injury prevention and post-injury stability for inversions, eversions and high ankle sprains
- Will not slip, slide or migrate
- #1 wrap in professional football & college sports



PowerFlex Catalog #	Size	Rolls/Case	Notes
3710XX-048	1" x 6 yds.	48	
3715XX-032	1.5" x 6 yds.	32	
3720XX-024	2" x 6 yds.	24	
3725XX-016	2.75" x 6 yds.	16	Available in
4725XX-048	2.75" x 6 yds.	48	non-latex
4730XX-048	3" x 6 yds.	48	(White)
3730XX-016	3" x 6 yds.	16	
3740XX-012	4" x 6 yds.	12	
3760XX-008	6" x 6 yds.	8	

Color Op	otions:								
		White (WH)	Yellow (YL)	Orange (OR)	Black (BK)	Teal (TE)	Tan (TN)	Red (RD)	Blue (BL)
N. Pink (NP)	Green (GN)	N. Green (NG)	Purple (PU)	Lt. Blue (LB)	Maroon (MR)	Grey (GY)	Navy (NV)	Black/Grey Camo(CK)	Woodland Camo(CW)

When ordering, please indicate catalog # and color code. For example, 3710WH-030 (XX = color code).

SYNTHETIC ADHESIVE SPORTS TAPE



POWERPRO

PowerPro® tape empowers stronger athletes by providing the most supportive adhesive taping system available. It is the first adhesive tape engineered with a non-stretch, moisture-stable substrate. Our most versatile adhesive tape, it pairs easily with any base.

From the first minute to final score, the PowerPro synthetic taping system sustains its integrity - even when wet.


- Stretches up to 37% less than cotton
- Sticks to anything PowerFlex, prewrap or skin
- Sustains its integrity when wet
- Hand-tears with ease and conforms around ankles
- Replaces traditional cloth athletic tape
- Provides support & stability during games



PowerPro Catalog #	Size	Rolls/Case
ACP165XX-1"	1" x 15 yds.	48
ACP165XX-1.5"	1.5" x 15 yds.	32
ACP165XX-2"	2" x 15 yds.	24

When ordering, please indicate catalog # and color code. For example, ACP165WH-1" (XX = color code).

STRETCH LAYER OPTIONS



POWERSPEED

PowerSpeed* is a high strength cohesive athletic wrap designed for speed taping or a stronger spat. It delivers maximum support for a modified taping system that won't slip/migrate from sweat or water. With easy hand-tear, PowerSpeed has quick application.

- 30% stronger than PowerFlex
- Supports sprains & strains



PowerFAST® is a strong, flexible adhesive stretch tape made with synthetic fabric. Unlike other adhesive stretch tapes, it won't absorb water so it remains thin & comfortable, rather than getting bulky. PowerFast also retains its elasticity for lasting support, plus its strong adhesive won't pull hair.

- Synthetic construction won't absorb water
- Thin and comfortable no bulk





PowerSpeed Catalog #	Size	Rolls/Case	PowerFAST Catalog #	Size	Rolls/Case
PS3720XX-024	2" x 6 yds.	24	4210XX-048	1" x 5 yds.	48
PS3725XX-016	2.75" x 6 yds.	16	4215XX-032	1.5" x 5 yds.	32
PS3730XX-016	3" x 6 yds.	16	4220XX-024	2" x 5 yds.	24
			4230XX-016	3" x 5 yds.	16
			4310XX-048	1" x 7.5 yds.	48
			4315XX-032	1.5" x 7.5 yds.	32
			4320XX-024	2" x 7.5 yds.	24
			4330XX-016	3" x 7.5 yds.	16

Color Options:

White Black

When ordering, please indicate catalog # and color code. For example, PS3720WH-024 (XX = color code).



POWERTAPE

PowerTape® is a cohesive athletic tape that sticks only to itself or to PowerFlex, creating one of the strongest combination taping systems in the field.

- 41 lbs. tensile strength
- Less than 1% stretch



VictoryTape® is an easier-to-tear and economical version of PowerTape to apply over PowerFlex for a long-lasting taping system.

- 38 lbs. tensile strength
- Sweat & water resistant

MOLESKIN

Moleskin is a tan-colored flannel with a napped cotton backing.







PowerTape Catalog #	Size	Rolls/Case	VictoryTape Catalog #	Size	Rolls/Case
ACP130XX-1"	1" x 15 yds.	48	ACP120XX-1"	1" x 15 yds.	48
ACP130XX-1.5"	1.5" x 15 yds.	32	ACP120XX-1.5"	1.5" x 15 yds.	32
ACP130XX-2"	2" x 15 yds.	24	ACP120XX-2"	2" x 15 yds.	24

Color O	ptions:
White	Black
 (WH)	(BK)

When ordering, please indicate catalog # and color code. For example, ACP130WH-1" (XX = color code).


Moleskin Catalog #	Size	Rolls/Case	Colors
041-02-024	2" x 25 yds.	24	tan only
041-03-016	3" x 25 yds.	16	tan only

DRESSINGS



POWERFLEX AFD

PowerFlex AFD® is an all-in-one foam pad dressing and cohesive bandage that helps get players back in the game. It provides fast, effective treatment for wounds, lacerations and turf burns.



- Individually packaged
- Pefect for game time kits



PowerFlex AFD Catalog #	Size	Rolls/Case
3722XX-016	2" x 2.5 yds.	16
3744XX-008	4" x 2.5 yds.	8

Color O	ptions:	
White (WH)	Black (BK)	
, ,	, ,	aso indic
catalog #	and color	r code.
Г	nla 27221	A/LI 016
	White (WH) When or catalog #	

APPLICATION AND REMOVAL TIPS

STEP ONE:

- Apply PowerFlex directly to the skin with proper unwind tension
- Remove all stretch when applying



STEP TWO:

 Cover PowerFlex with PowerPro, PowerTape, or VictoryTape



POST-INJURY TURF TOE/ANKLE COMBINATION



BASIC PREVENTIVE WRIST/HAND/THUMB



POST-INJURY TURF TOE



REGULAR OR BASIC PREVENTIVE ANKLE



POST-INJURY ELBOW



SPEED SPATTING



Milliken

Additional resources for our sports medicine options including product information sheets, application videos and case studies are available at healthcare.milliken.com.





