

Sharppoint[™] **PLUS**
Surgical Sutures

SUTURE CATALOG

Choose Excellence, Expect Value



**SURGICAL
SPECIALTIES**
CORPORATION[™]

Choose Excellence, Expect Value

Surgical Specialties Corporation is the world's largest independent suture manufacturer. While headquartered in Westwood, Massachusetts, SSC is a global company with development and manufacturing facilities in the United States, Germany, Mexico, China, and the United Kingdom.

We design and manufacture high performance wound closure products including QUILL™, Sharpoint™, and LOOK™ brand sutures.

We are pleased to offer the new Sharpoint™ PLUS suture portfolio featuring German engineered needles.

Our philosophy is that performance doesn't have to be expensive.



Table of Contents

Suture Material Selection	3-4
Advanced Needle Technology	5
Needle Anatomy and Terminology	6
Needle Points and Cutting Edges	7
Needle Charts.....	8-10
 Absorbable Sutures	
Plain Gut.....	11
Chromic Gut	12-13
PolySyn FA™ (polyglycolic acid)	14
Monoderm™ (glycolide and ε-caprolactone)	15-16
PolySyn™ (polyglycolic acid).....	17-21
PDO (polydioxanone)	22
 Non-Absorbable Sutures	
Silk.....	23-25
Nylon™	26-28
Polypropylene.....	29-32
Polyviolen™ (polyester) - Uncoated	33
Polyviolen™ XL (polyester)- Coated.....	34
 Cross Reference Guide	
Ethicon Cross Reference.....	35-41
Covidien Cross Reference.....	42-45






SharpPoint **PLUS**
Surgical Sutures™

Guide to Suture Materials

Absorbable


Material	Material Description	BSR Profile	Mass Absorption
GUT  Comparable to: Surgical Gut Suture - Plain	Twisted • Plain & Chromic • Packaged Dry or in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
CHROMIC GUT  Comparable to : Surgical Gut Suture - Chromic	Twisted • Plain & Chromic • Packaged Dry or in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
POLYSYN FA™ (polyglycolic acid)  Comparable to: VICRYL RAPIDE™	Braided • Undyed (clear) • Coated Common uses: General Closure, OB/Gyn, Urology	In vivo strength retention: 48% at 7 days 0% at 14 days	Essentially complete between 49 and 63 days
MONODERM™ (polyglycolic acid and ε-caprolactone)  Comparable to: MONOCRYL™	Monofilament • Undyed (clear) • Dyed (violet) Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 42% to 76% at 7 days 36% to 52% at 14 days	Essentially complete by 90 days
POLYSYN™ (polyglycolic acid)  Comparable to: VICRYL™	Braided • Undyed (clear) • Dyed (violet) Monofilament (9-0 & 10-0 only) • Undyed (clear) • Dyed (violet) • Coated & Uncoated Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 50% at 14 days 20% at 21 days	Essentially complete between 50 and 90 days
PDO (polydioxanone)  Comparable to: PDS II™	Monofilament • Dyed (violet) Common uses: Plastic Surgery, Orthopaedic, OB/Gyn, Gastrointestinal	In vivo strength retention: 80% to 90% at 14 days 60% to 82% at 28 days 47% to 79% at 42 days	Essentially complete between 180 and 220 days

Non-Absorbable

Material	Material Description	Qualities
SILK  Comparable to: PERMA-HAND™ Silk	Braided • Undyed (white) • Dyed (black) • Coated Common uses: General & Skin Closure, Gastrointestinal, Cardiovascular Surgery, Plastic Surgery, Ophthalmology	• Excellent handling and tying characteristics
NYLON™  Comparable to: ETHILON™	Monofilament • Undyed (clear) • Dyed (blue and black) • Packaged dry or in wetting solution Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery, Gastrointestinal	• Passes easily through tissue
POLYPROPYLENE  Comparable to: PROLENE™	Monofilament • Dyed (blue) Common uses: General & Skin Closure, Neurosurgery, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery	• Inert • Easy to handle • Passes easily through tissue
POLYVIOLENE™ (UNCOATED) (poly (ethylene-terephthalate))  Comparable to: MERSILENE™	Braided • Undyed (white) • Dyed (green) • Uncoated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security
POLYVIOLENE™ XL (COATED)  Comparable to: ETHIBOND™	Braided • Undyed (white) • Dyed (green) • Coated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security

Advanced Needle Technology

SharpTM needle performance is an output of its design. Materials, geometry, and coatings work together to produce a high performance surgical instrument.




Proprietary Stainless Steel Alloy


SharpTM needles are forged from an alloy of Iron, Nickel, and Chromium. Each element serves a role in delivering strength, ductility, and hardness.



Strength




Ductility





Hardness

Robotic Manufacturing

SharpTM needles are manufactured using German designed robotics to ensure quality and repeatability.



Proven Performance



*Testing performed by independent lab.
Acculine Engineering; Waterbury, CT*

Needle Terminology

Length: Distance measured along the circumference of the needle, from the point to the swage. It is depicted on the packaging.

Shape: Different surgeries or surgical techniques demand different needle shapes.

Point: The sharpened part of the needle that first penetrates the tissue. The point consists of the area from the sharpened tip of the needle to the widest cross section of the needle body. The point determines ease of penetration and the initial size and shape of the hole.


Swage: Connection point of the suture and the needle. Micro teeth grip the suture material and prevent it from slipping out.

Wire Diameter: Gauge or thickness of the wire originally used to form the needle body.


Body: The part between the point and the swage of the needle.

Included Angle (Curvature): The fraction of the circle included in the curve of the needle. This is depicted on the packaging.


Needle Anatomy




Needle Shapes and General Uses




Straight (aka: Keith): Skin Closure




Bi-curve: Vasectomy Reversal, Ophthalmic




1/4 Circle: Ophthalmic



1/8 Circle: Ophthalmic



3/8 Circle: Skin Closure, Ophthalmic



1/2 Circle: Skin Closure, GI, Cardiovascular, Ophthalmic

Needle Points And Cutting Edges

Needle Type	Point Geometry	Description	Commonly Used For
Conventional Cutting		Triangular point with cutting edge on the inner curvature	General skin closure, subcutaneous tissue, sometimes for ophthalmic surgery, plastic or reconstructive surgery
Precision Conventional Cutting		Triangular point with cutting edge on the inner curvature and precision-honed for added sharpness	General skin closure, subcutaneous tissue, sometimes for ophthalmic surgery, plastic or reconstructive surgery
Reverse Cutting		Triangular point with cutting edge on the outer curvature	Skin closure, ligamentous or fibrous tissues
Precision Reverse Cutting		Triangular point with cutting edge on the outer curvature and precision-honed for added sharpness	Primarily used in plastic surgery for delicate procedures and where a good cosmetic result is desired
Taper Point		Round shaft, straight or curved, taper point with no cutting edge	Soft tissue closure such as fascia, vascular, gastrointestinal, and most soft tissue below the skin
Taper Cutting		Four (4) cutting edges with a taper shaft	For penetration through tough tissue
Spatula		Spatula needle with the cutting edge on the inner or outer curvature	Ophthalmic and microsurgery
Center Point Spatula		Two cutting edges in a horizontal plane	Ophthalmic surgery for muscle and retinal repair. Also used for delicate plastic surgery
Straight Cutting		Reverse cutting tip on a straight shaft	General skin closure
M.E.T.™ Point		Triangular faceted point tapering to a round body	Microsurgery
UltraGlide™		Six (6) faceted sides and two (2) cutting edges	Ophthalmic surgery for cornea and sclera closure
Vas Cutting		Cutting edges aid in penetration of multiple layers of tissue	Microsurgery and urology

Center Point Spatula Needles



Surgical Specialties	Ethicon™	Covidien™	
1/4 Circle			
DO5	RD-1	SS-24	8mm

M.E.T.™ Point Needles



Surgical Specialties	Ethicon™	Covidien™	
Bi-curve			
BRM5	•	•	5mm
1/4 Circle			
VRM4	V75-3	•	4mm
VRM5	•	•	5mm
VRM6	•	•	6mm
3/8 Circle			
DRM4	•	MVK-70-3	4mm
DRM5	•	MVK-100-4	5mm
DRM6	•	•	6mm
1/2 Circle			
HRM4	•	•	4mm

Precision Conventional Cutting Needles



Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
PC-1	PC-1	P-13	13mm
PC-3	PC-3	PC-11	16mm
PC-5	PC-5	P-12	19mm



Needle Chart

Precision Reverse Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
G-6	G-6	HE-6	7mm
P-1	P-1	P-10	11mm
G-1	G-1	G-1	11mm
P-3	P-3	P-13	13mm
G-3	G-3	G-3	13mm
PS-3	PS-3	P-11	16mm
PS-2	PS-2	P-12	19mm
PS-1	PS-1	P-14	24mm
1/2 Circle			
G-7	G-7	HE-7	8mm

Reverse Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
C-2	C-2	C-1	11mm
C-2	C-2	C-1	12mm
C-3	C-3	C-12	14mm
FS-3	FS-3	C-12	16mm
FS-2	FS-2	C-13	18mm
FS-1	FS-1	C-14	24mm
FS	FS	C-15	26mm
FSL	FSL	C-16	30mm
FSLX	FSLX	C-17	36mm
FSLX	FSLX	C-17	40mm
DS90	•	•	90mm
1/2 Circle			
M-1	M-1	C-21	15mm
J-1	J-1	C-22	20mm
X-1	X-1	C-23	24mm
CP-2	CP-2	GS-10	26mm
CP-1	CP-1	GS-11	37mm
CP	CP	GS-12	40mm

Spatula Needles

Surgical Specialties	Ethicon™	Covidien™	
Bi-curve			
BSL5	TG6-C	•	5.51mm
1/8 Circle			
SSL5	•	•	5.51mm
SSL15	•	•	14.99mm
1/4 Circle			
VSL7	TG100-8	SE100-8	6.6mm
VSP8	S-29, S-24	SD-1, SS-24	8mm
VCL8	S-14	SS-14	8.51mm
3/8 Circle			
DSL6	TG140-6	SE140-6	6.15mm
DSL7	TG140-8	SE140-8	6.6mm
1/2 Circle			
HSL5	TG160-6	SE160-6	5.51mm
HSL6	TG160-8	SE160-8	5.54mm
HSL7	•	•	7mm


Taper Cutting Needles



Surgical Specialties Ethicon™ Covidien™


3/8 Circle

DRM4 BV-100-3 MVK-70-3
 DRM5 BV-100-4 •
 DRM6 BV-130-5 MV-135-5
 • CC-1 •
 • RB-1 CV-23



1/2 Circle

• V-5 KV-5
 HRT17 V-5 KV-5
 HRT31 V-34 KV-34
 HRT36 V-34 KV-34
 HRT38 V-34 KV-34



Taper Point Needles



Surgical Specialties Ethicon™ Covidien™


1/4 Circle

VR5 V75-3 •




3/8 Circle

DR4 BV75-3 MV70-3
 DR5 BV130-4 MV135-4
 DR6 BV130-5 MV135-5
 DR12 C-1 CV-1
 DR20 SH-1 CV-24



1/2 Circle

HR13 TF CV-22
 HR18 RB-1 CV-23
 HR22 SH-1 CV-25
 HR26 SH V-20
 • CT-2 GS-22
 HR30 • V-30
 HR36 MH V-26
 HRX36 CT-1 GS-21
 HRX40 CT GS-24X



Straight Cutting Needles



Surgical Specialties Ethicon™ Covidien™

GS60

KS

SC-2

60mm

UltraGlide™ Needles



Surgical Specialties Ethicon™ Covidien™

Bi-curve

BGL5

CSB-6

•

4.88mm

1/4 Circle

VGL7

CS90-6

•

7mm

3/8 Circle

DGL6

CS140-6

SE-140-6

6.15mm

DGL7

CS140-6

SE-140-6

6.60mm

1/2 Circle

HGL5

CS160-6

SE-160-6

5.51mm

HGL7

CS175-6

•

6.99mm

Vas Cutting Needles



Surgical Specialties Ethicon™ Covidien™

1/2 Circle

HSV6

•

•

6mm

Plain Gut

Precision Reverse Cutting				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
G-1 11mm 3/8 Circle		Plain Gut	18"/45cm	B774N					
G-1 11mm 3/8 Circle		Plain Gut	18"/45cm	B770N					
P-3 13mm 3/8 Circle		Plain Gut	18"/45cm		B686N	B1644N			
G-3 13mm 3/8 Circle		Plain Gut	18"/45cm		B772N				
PS-2 19mm 3/8 Circle		Plain Gut	18"/45cm			B1627N			

Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle		Plain Gut	27"/70cm			B821N	B822N		

Straight Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
SC-1 13mm Straight 		Plain Gut	18"/45cm			B1824N			
SC-1 13mm Straight 		Plain Gut	18"/45cm			B1828N			

N=box of 12; D=box of 36

Chromic Gut

Precision Reverse Cutting		▽		Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
G-1 11mm 3/8 Circle		Chromic Gut	18"/45cm	C796N					
G-1 11mm 3/8 Circle		Chromic Gut	18"/45cm	C790N					
P-3 13mm 3/8 Circle		Chromic Gut	18"/45cm		C687N	C1654N			
G-3 13mm 3/8 Circle		Chromic Gut	18"/45cm		C792N	C793N			
PS-3 16mm 3/8 Circle		Chromic Gut	18"/45cm		C1636N				
PS-4 16mm 1/2 Circle		Chromic Gut	18"/45cm			C1643N			
PS-2 19mm 3/8 Circle		Chromic Gut	18"/45cm			C1637N			
			27"/70cm				C1638N		

Reverse Cutting		▽		Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle		Chromic Gut	27"/70cm			C635N	C636N		

N=box of 12; D=box of 36

Chromic Gut

Taper Point 		Size							
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
T28 18mm 1/2 Circle RB-1 CV-23		Chromic Gut	14"/35cm		C0590N				
HR22/T4 22mm 1/2 Circle SH-1 CV-25		Chromic Gut	14"/35cm			C181N			
			30"/75cm				C182N		
HR26 26mm 1/2 Circle SH V-20		Chromic Gut	30"/75cm			C121N	C122N	C123N	
HR30 30mm 1/2 Circle V-20		Chromic Gut	30"/75cm					C0244N	
HR36 36mm 1/2 Circle MH V-26		Chromic Gut	30"/75cm						C128N

N=box of 12; D=box of 36

Precision Reverse Cutting

Size

		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3		Undyed Braided	18"/45cm		M493N	M494N			
13mm 3/8 Circle									
PS-2		Undyed Braided	18"/45cm			M496N	M497N		
19mm 3/8 Circle			27"/70cm			M426N			

Taper Cutting

Size

		Material	Length	6-0	5-0	4-0	3-0	2-0	0
HRT36		Undyed Braided	36"/90cm					M917N	
36mm 1/2 Circle V-34 KV-34									

Taper Point

Size

		Material	Length	6-0	5-0	4-0	3-0	2-0	0
HR17		Undyed Braided	27"/70cm			M214N			
17mm 1/2 Circle RB-1 CV-23									

N=box of 12; D=box of 36

Precision Conventional Cutting

Precision Conventional Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1 13mm 3/8 Circle		Undyed Monofilament	18"/45cm		Y834N				
PC-3 16mm 3/8 Circle		Undyed Monofilament	18"/45cm		Y844N	Y845N			
PC-5 19mm 3/8 Circle		Undyed Monofilament	18"/45cm			Y823N			

Precision Reverse Cutting

Precision Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3  13mm 3/8 Circle	Undyed Monofilament	18"/45cm	Y492N	Y493N	Y494N				
	Violet Monofilament	18"/45cm		Y463N	Y464N				
PS-2  19mm 3/8 Circle	Undyed Monofilament	18"/45cm		Y495N	Y496N	Y497N			
		27"/70cm			Y426N	Y427N			
	Violet Monofilament	18"/45cm			Y513N				
PS-1  24mm 3/8 Circle	Undyed Monofilament	18"/45cm			Y682N				
		27"/70cm			Y935N	Y936N			

N=box of 12; D=box of 36

Reverse Cutting ▽

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
CP-1  36mm 1/2 Circle	Undyed Monofilament	27"/70cm					Y266N	

Taper Point ⊙

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
HR13  13mm 1/2 Circle TF CV-22	Undyed Monofilament	27"/70cm	Y432N	Y433N	Y434N			
HR26  26mm 1/2 Circle CT-2 NCT-2	Undyed Monofilament	27"/70cm					Y417D	

PolySyn™/Polyglycolic Acid

Precision Conventional Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-3 16mm 3/8 Circle		Undyed Braided	18"/45cm			G845N			
PC-5 19mm 3/8 Circle		Undyed Braided	18"/45cm				G824N		

Precision Reverse Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-1 11mm 3/8 Circle		Undyed Braided	18"/45cm	G489N					
P-3 13mm 3/8 Circle		Undyed Braided	18"/45cm	G492N	G493N	G494N			
		Violet Braided	18"/45cm		G463N	G464N			
PS-3 16mm 3/8 Circle		Undyed Braided	18"/45cm		G500N				
PS-4 16mm 1/2 Circle		Undyed Braided	18"/45cm			G507N			
PS-2 19mm 3/8 Circle		Undyed Braided	18"/45cm		G495N	G496N	G497N		
			27"/70cm			G426N	G427N	G428N	
PS-1 24mm 3/8 Circle		Undyed Braided	18"/45cm				G683N		
			27"/70cm			G935N	G936N		

N=box of 12; D=box of 36

PolySyn™/Polyglycolic Acid

Reverse Cutting ▼			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
FS-2  19mm 3/8 Circle	Undyed Braided	27"/70cm		G422N	G423N			
	Violet Braided	18"/45cm		G392N	G393N			
FS-1  24mm 3/8 Circle	Undyed Braided	27"/70cm		G441N	G442N	G443N		
CP-2  26mm 1/2 Circle	Undyed Braided	27"/70cm				G869N		
CP-1  36mm 1/2 Circle	Undyed Braided	27"/70cm				G266N	G267N	G268N
OS-6  36mm 1/2 Circle	Undyed Braided	27"/70cm						G535N

N=box of 12; D=box of 36

PolySyn™/Polyglycolic Acid

Spatula



			Size			
	Material	Length	6-0	5-0	4-0	3-0
VSL7  6.60mm 1/4 Circle 6.60mm TG100-8 SE-100-8	Violet Braided	18"/45cm	G544N			
VCL8  8.51mm 1/4 Circle SG-14 SS-14	Undyed Braided	18"/45cm	G670N	G671N		
	Violet Braided	8"/20cm	G590N			
VSP8  7.62mm 1/4 Circle S-29 SD-1	Violet Braided	12"/30cm	G556N			
		18"/45cm	G555N			
VSP08  8mm 1/4 Circle S-24 SS-24	Violet Braided	12"/30cm	G552N	G553N		

Taper Cutting



			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
V-34  36mm 1/2 Circle	Undyed Braided	36"/90cm					G917N	

N=box of 12; D=box of 36

Taper Point 		Size							
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
HR13 13mm 1/2 Circle TF CV-22		Undyed Braided	27"/70cm		G433N				
		Violet Braided	18"/45cm	G981N					
HR17 17mm 1/2 Circle RB-1 CV-23		Undyed Braided	27"/70cm			G714N			
		Violet Braided	30"/75cm			G0304N			
T14 22mm 1/2 Circle SH-1 CV-24		Violet Braided	27"/70cm			G310N	G311N		
HR26 26mm 1/2 Circle SH V-20		Undyed Braided	27"/70cm				G416N		
		Violet Braided	27"/70cm				G316N	G317N	
T5 26.5mm 1/2 Circle CT-2 GS-22		Undyed Braided	27"/70cm				G232N	G269N	
		Violet Braided	27"/70cm				G332N	G333N	
HR30 30mm 1/2 Circle • V-30		Violet Braided	30"/75cm				G0322N	G323N	

N=box of 12; D=box of 36

PolySyn™/Polyglycolic Acid

Taper Point 			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
T8  36.6mm 1/2 Circle CT-1 GS-21	Undyed Braided	18"/45cm				G839N		
		27"/70cm					G260N	
		36"/90cm				G945N		
	Violet Braided	27"/70cm				G339N	G340N	G341N
		36"/90cm				G354N	G0518N	
T10  40mm 1/2 Circle CT GS-24	Violet Braided	27"/70cm					G352N	G353N
		36"/90cm					G358N	G359N

N=box of 12; D=box of 36

Precision Reverse Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3 13mm 3/8 Circle 		Violet Monofilament	18"/45cm		R493N	R494N			
PS-2 19mm 3/8 Circle 		Violet Monofilament	18"/45cm			R513N	R497N		

Reverse Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle 		Violet Monofilament	27"/70cm			R422N	R423N		
FS-1 24mm 3/8 Circle 		Violet Monofilament	27"/70cm					R443N	

Taper Point

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
HR17 17mm 1/2 Circle RB-1 NRB-1 		Violet Monofilament	27"/70cm				R305D		
T5 26mm 1/2 Circle SH/CT-2 NSH/NC-2 		Violet Monofilament	27"/70cm				R316D	R317D/ R333D	R334D
HR3 36mm 1/2 Circle CT-1 NCT-1 		Violet Monofilament	27"/70cm					R339D	R340D

N=box of 12; D=box of 36

Precision Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-1 11mm 3/8 Circle 		Black Braided	18"/45cm	D639N					
G-1 11mm 3/8 Circle 		Black Braided	18"/45cm	D780N					
P-3 13mm 3/8 Circle 		Black Braided	18"/45cm		D640N	D641N			
G-3 13mm 3/8 Circle 		Black Braided	18"/45cm			D783N			
PS-1 24mm 3/8 Circle 		Black Braided	18"/45cm				D1684N		
PS 26mm 3/8 Circle 		Black Braided	18"/45cm					D1588N	
PSL 30mm 3/8 Circle 		Black Braided	18"/45cm						D580N
			30"/75cm						D590N

Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-1 24mm 3/8 Circle 		Black Braided	18"/45cm				D684N		
FS 26mm 3/8 Circle 		Black Braided	18"/45cm					D685N	

N=box of 12; D=box of 36

Reverse Cutting



			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
FSL  30mm 3/8 Circle	Black Braided	18"/45cm					D678N	
		30"/75cm					D680N	

Spatula



			Size		
	Material	Length	8-0	7-0	6-0
VSP8  8mm 1/4 Circle TG100-8 SE-100-8	Black Braided	18"/45cm			D1780N
DSL7  6.15mm 3/8 Circle TG140-8 SE-140-8	Black Braided	12"/30cm	D2526N		
		18"/45cm		D7733N	

Straight Cutting



			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
GS60  60mm Straight KS SC-2	Black Braided	30"/75cm					D623N	

N=box of 12; D=box of 36

Silk

Taper Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
DRT14  13mm 3/8 Circle CC-1 KZ-11	Black Braided	18"/45cm	D7729N					

Taper Point

			Size						
	Material	Length	6-0	5-0	4-0	3-0	2-0	0	1
DR12  13mm 3/8 Circle C-1 CV-1	Black Braided	18"/45cm	D707N		D7734N				
HR26  26mm 1/2 Circle SH V-20	Black Braided	30"/75cm				D832N			
HR36  36mm 1/2 Circle MH V-26	White Braided	30"/75cm							D2543N

Pre-Cut Sutures

		Size				
	Length	0	1	2	3	4
Silk, Black Braided Strands	8 - 24"/60cm					D2578N
	13 - 24"/60cm	D76N		D0492N		

N=box of 12; D=box of 36

Precision Conventional Cutting

Precision Conventional Cutting 				Size					
Needle		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1  13mm 3/8 Circle		Black Monofilament	18"/45cm	A1856N	A1855N				
		Black Monofilament Pliable	18"/45cm		A1955N				
PC-3  16mm 3/8 Circle		Black Monofilament	18"/45cm		A1865N	A1864N			
		Black Monofilament Pliable	18"/45cm		A1965N				
PC-5  19mm 3/8 Circle		Black Monofilament	18"/45cm			A1894N	A1893N		
		Black Monofilament Pliable	18"/45cm			A1994N			

Precision Reverse Cutting

			Size					
Needle		Material	Length	6-0	5-0	4-0	3-0	2-0
P-1 11mm 3/8 Circle 		Black Monofilament	18"/45cm	A697N				
P-3 13mm 3/8 Circle 		Black Monofilament	18"/45cm	A1698N	A698N	A699N		
		Undyed Monofilament	18"/45cm		A690N	A691N		
PS-4 16mm 1/2 Circle 		Black Monofilament	18"/45cm			A1662N		
PS-3 16mm 3/8 Circle 		Black Monofilament	18"/45cm	A1665N	A1668N			

N=box of 12; D=box of 36

Nylon

Precision Reverse Cutting



Size

		Material	Length	6-0	5-0	4-0	3-0	2-0
PS-2 19mm 3/8 Circle		Black Monofilament	18"/45cm		A1666N	A1667N	A1669N	A593N
		Blue Monofilament	18"/45cm			AG667N		
PS-1 24mm 3/8 Circle		Black Monofilament	18"/45cm				A1663N	
PS 26mm 3/8 Circle		Black Monofilament	18"/45cm					A585N
PSL 30mm 3/8 Circle		Blue Monofilament	30"/75cm				A1691N	
PSLX 36mm 3/8 Circle		Black Monofilament	30"/75cm				A1683N	A1697N

Reverse Cutting



Size

		Material	Length	5-0	4-0	3-0	2-0	0	1
FS-2 19mm 3/8 Circle		Black Monofilament	18"/45cm	A661N	A662N				
FS-1 24mm 3/8 Circle		Black Monofilament	18"/45cm		A1629N	A663N			
			30"/75cm			A669N			
FS 26mm 3/8 Circle		Black Monofilament	18"/45cm				A664N		

N=box of 12; D=box of 36

Reverse Cutting ▽			Size					
Needle	Material	Length	5-0	4-0	3-0	2-0	0	1
FSL  30mm 3/8 Circle	Black Monofilament	30"/75cm			A1671N			
FSLX  36mm 3/8 Circle	Black Monofilament	30"/75cm			A1673N	A1674N		

Straight Cutting ▽			Size				
Needle	Material	Length	6-0	5-0	4-0	3-0	2-0
SC2  60mm Straight KS SC-2	Black Monofilament	30"/75cm				A627N	

Polypropylene

Precision Conventional Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1 13mm 3/8 Circle		Blue Monofilament	18"/45cm	J8617N	J8618N				
PC-3 16mm 3/8 Circle		Blue Monofilament	18"/45cm		J8635N	J8634N			
PC-5 19mm 3/8 Circle		Blue Monofilament	18"/45cm			J8631N	J8632N		

Precision Reverse Cutting

				Size				
		Material	Length	7-0	6-0	5-0	4-0	3-0
P-1 11mm 3/8 Circle		Blue Monofilament	18"/45cm		J8697N			
P-3 13mm 3/8 Circle		Blue Monofilament	18"/45cm		J8695N	J8698N	J8699N	
PS-3 16mm 3/8 Circle		Blue Monofilament	18"/45cm		J8680N	J8681N		
PS-2 19mm 3/8 Circle		Blue Monofilament	18"/45cm			J8686N	J8682N	J8687N
PS-1 24mm 3/8 Circle		Blue Monofilament	18"/45cm					J8663N

Straight Cutting

				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
GS60 60mm Straight KS SC-2		Blue Monofilament	30"/75cm					J8623N	

N=box of 12; D=box of 36

Reverse Cutting


				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle		Blue Monofilament	18"/45cm		J8661N	J8683N	J8665N		
FS-1 24mm 3/8 Circle		Blue Monofilament	18"/45cm				J8684N		
FS 26mm 3/8 Circle		Blue Monofilament	18"/45cm					J8685N	

Taper Point

				Size				
		Material	Length	7-0	6-0	5-0	4-0	3-0
BV-1 9.3mm 3/8 Circle		Blue Monofilament	18"/45cm	J8701D	J8806D			
			24"/60cm	J8702D	J8805D			
			30"/75cm	J8703D	J8709D			
C-1 13mm 3/8 Circle		Blue Monofilament	18"/45cm		J8718D	J8717D		
			24"/60cm		J8726D	J8725D		
			30"/75cm		J8706D			
			36"/90cm			J8720D		
RB-2 13mm 1/2 Circle		Blue Monofilament	18"/45cm		J8714D	J8713D		
			24"/60cm			J8554D		
			30"/75cm		J8711D	J8710D		

N=box of 12; D=box of 36

Polypropylene

Taper Point 				Size				
		Material	Length	6-0	5-0	4-0	3-0	2-0
BB 17mm 3/8 Circle		Blue Monofilament	36"/90cm		J8580D	J8581D		
RB-1 17mm 1/2 Circle		Blue Monofilament	18"/45cm		J8756D			
			24"/60cm		J8555D			
			36"/90cm		J8556D	J8557D	J8558D	
SH-1 22mm 1/2 Circle		Blue Monofilament	30"/75cm			J8526D	J8762D	
SH 26mm 1/2 Circle		Blue Monofilament	36"/90cm			J8521D	J8522D	J8523D
			48"/122cm				J8534D	J8533D
MH 36mm 1/2 Circle		Blue Monofilament	36"/90cm				J8842D	J8843D
			48"/122cm					J8853D

N=box of 12; D=box of 36

Taper Cutting

Taper Cutting 				Size					
		Material	Length	7-0	6-0	5-0	4-0	3-0	2-0
CC 9.3mm 3/8 Circle 		Blue Monofilament	24"/60cm	J8704D	J8807D				
CC-1 13mm 3/8 Circle 		Blue Monofilament	30"/75cm		J8707D				
V-5 17mm 1/2 Circle 		Blue Monofilament	36"/90cm				J8935D		
V-7 26mm 1/2 Circle 		Blue Monofilament	36"/90cm				J8975D	J8976D	J8977D

Polyviolene™ / Polyester - Uncoated

Reverse Cutting ▽

			Size					
	Material	Length	7-0	6-0	5-0	4-0	3-0	2-0
C22  16mm 3/8 Circle FS-3 C-12	Green Braided Uncoated	18"/45cm			F681N			
DS18  18mm 3/8 Circle FS-2 C-13	Green Braided Uncoated	18"/45cm				F683N		
	White Braided Uncoated	18"/45cm				F633N		
DS24  24mm 3/8 Circle FS-1 C-14	White Braided Uncoated	18"/45cm					F647D	
DS26  26mm 3/8 Circle FS C-15	White Braided Uncoated	18"/45cm						F665D

Side Cutting Lancet ▽

			Size					
	Material	Length	7-0	6-0	5-0	4-0	3-0	2-0
VSP8  8mm 1/4 Circle S-24 SS-24	White Braided Uncoated	18"/45cm			F893N			

Taper Point ●

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
HRT17  17mm 1/2 Circle V-5 KV-5	Green Braided Uncoated	12"/30cm Loop		F951N	F2600N	F872N		

N=box of 12; D=box of 36

Polyviolene™ XL / Polyester - Coated

Reverse Cutting ▽

				Size				
	Material	Length	5-0	4-0	3-0	2-0	0	2
OS-4  22mm 1/2 Circle	Green Braided	30"/75cm						E519N

Taper Point ◎

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
RB-1  17mm 1/2 Circle	Green Braided	30"/75cm			EX551D			
		36"/90cm		EX556D	EX557D	EX558D		
SH-1  22mm 1/2 Circle	Green Braided	36"/90cm				EX762D	EX763D	
SH  26mm 1/2 Circle	Green Braided	30"/75cm					EX563D	
		36"/90cm					EX523D	
	White Braided	36"/90cm					EX513D	

Taper Cutting ◎

			Size					
Needle	Material	Length	6-0	5-0	4-0	3-0	2-0	0
V-4  17mm 3/8 Circle	Green Braided	30"/75cm					EX927D	
V-7  26mm 1/2 Circle	Green Braided	30"/75cm					EX977D	

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
1588H	D1588N	Silk 2-0 Black 18" PS Fine Reverse Cutting 3/8C 26mm
1627H	B1627N	Plain Gut 4-0 18" PS-2 Fine Reverse Cutting 3/8C 19mm
1629H	A1629N	Nylon 4-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
1636G	C1636N	Chromic Gut 5-0 18" PS-3 Fine Reverse Cutting 3/8C 16mm
1637G	C1637N	Chromic Gut 4-0 18" PS-2 Fine Reverse Cutting 3/8C 19mm
1638H	C1638N	Chromic Gut 3-0 27" PS-2 Fine Reverse Cutting 3/8C 19mm
1643G	C1643N	Chromic Gut 4-0 18" PS-4 Fine Reverse Cutting 1/2C 16mm
1644G	B1644N	Plain Gut 4-0 18" P-3 Fine Reverse Cutting 3/8C 13mm
1654G	C1654N	Chromic Gut 4-0 18" P-3 Fine Reverse Cutting 3/8C 13mm
1662G	A1662N	Nylon 4-0 Black 18" PS-4 Fine Reverse Cutting 1/2C 16mm
1663G	A1663N	Nylon 3-0 Black 18" PS-1 Fine Reverse Cutting 3/8C 24mm
1665G	A1665N	Nylon 6-0 Black 18" PS-3 Fine Reverse Cutting 3/8C 16mm
1666G	A1666N	Nylon 5-0 Black 18" PS-2 Fine Reverse Cutting 3/8C 19mm
1667G	A1667N	Nylon 4-0 Black 18" PS-2 Fine Reverse Cutting 3/8C 19mm
1668G	A1668N	Nylon 5-0 Black 18" PS-3 Fine Reverse Cutting 3/8C 16mm
1669H	A1669N	Nylon 3-0 Black 18" PS-2 Fine Reverse Cutting 3/8C 19mm
1671H	A1671N	Nylon 3-0 Black 30" FSL Reverse Cutting 3/8C 30mm
1673H	A1673N	Nylon 3-0 Black 30" FSLX Reverse Cutting 3/8C 36mm
1674H	A1674N	Nylon 2-0 Black 30" FSLX Reverse Cutting 3/8C 36mm
1683H	A1683N	Nylon 3-0 Black 30" PSLX Reverse Cutting 3/8C 36mm
1684G	D1684N	Silk 3-0 Black 18" PS-1 Fine Reverse Cutting 3/8C 24mm
1691H	A1691N	Nylon 3-0 Black 30" PSL Reverse Cutting 3/8C 30mm
1697H	A1697N	Nylon 2-0 Black 30" PSLX Reverse Cutting 3/8C 36mm
1698G	A1698N	Nylon 6-0 Black 18" P-3 Fine Reverse Cutting 3/8C 13mm
1824H	B1824N	Plain Gut 4-0 18" SC-1 Reverse Cutting Straight 13mm
1828H	B1828N	Plain Gut 4-0 18" SC-1 Reverse Cutting Straight 13mm Double Armed
1855G	A1855N	Nylon 5-0 Black 18" PC-1 Fine Conventional Cutting 3/8C 13mm
1856G	A1856N	Nylon 6-0 Black 18" PC-1 Fine Conventional Cutting 3/8C 13mm
1864G	A1864N	Nylon 4-0 Black 18" PC-3 Fine Conventional Cutting 3/8C 16mm
1865G	A1865N	Nylon 5-0 Black 18" PC-3 Fine Conventional Cutting 3/8C 16mm
1893G	A1893N	Nylon 3-0 Black 18" PC-5 Fine Conventional Cutting 3/8C 19mm
1894G	A1894N	Nylon 4-0 Black 18" PC-5 Fine Conventional Cutting 3/8C 19mm
1955G	A1955N	Nylon (pliable) 5-0 Black 18" PC-1 Fine Conventional Cutting 3/8C 13mm

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
1965G	A1965N	Nylon (pliable) 5-0 Black 18" PC-3 Fine Conventional Cutting 3/8C 16mm
1994G	A1994N	Nylon (Pliable) 4-0 Black 18" PC-5 Fine Conventional Cutting 3/8C 19mm
580H	D580N	Silk 0 Black 18" PSL Reverse Cutting 3/8C 30mm
585H	A585N	Nylon 2-0 Black 18" PS Fine Reverse Cutting 3/8C 26mm
590H	D590N	Silk 0 Black 30" PSL Reverse Cutting 3/8C 30mm
593H	A593N	Nylon 2-0 Black 18" PS-2 Fine Reverse Cutting 3/8C 19mm
623H	D623N	Silk, 2-0, Black, 30", GS60, Straight Straight Cutting 60mm
627H	A627N	Nylon, 3-0, Black, 30", GS60, Straight Straight Cutting 60mm
635H	C635N	Chromic Gut 4-0 27" FS-2 Reverse Cutting 3/8C 19mm
636H	C636N	Chromic Gut 3-0 27" FS-2 Reverse Cutting 3/8C 19mm
639G	D639N	Silk 6-0 Black 18" P-1 Fine Reverse Cutting 3/8C 11mm
640G	D640N	Silk 5-0 Black 18" P-3 Fine Reverse Cutting 3/8C 13mm
641G	D641N	Silk 4-0 Black 18" P-3 Fine Reverse Cutting 3/8C 13mm
661H	A661N	Nylon 5-0 Black 18" FS-2 Reverse Cutting 3/8C 19mm
662G	A662N	Nylon 4-0 Black 18" FS-2 Reverse Cutting 3/8C 19mm
663G	A663N	Nylon 3-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
664G	A664N	Nylon 2-0 Black 18" FS Reverse Cutting 3/8C 26mm
669H	A669N	Nylon 3-0 Black 30" FS-1 Reverse Cutting 3/8C 24mm
678G	D678N	Silk 0 Black 18" FSL Reverse Cutting 3/8C 30mm
680H	D680N	Silk 0 Black 30" FSL Reverse Cutting 3/8C 30mm
684G	D684N	Silk 3-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
685G	D685N	Silk 2-0 Black 18" FS Reverse Cutting 3/8C 26mm
686G	B686N	Plain Gut 5-0 18" P-3 Fine Reverse Cutting 3/8C 13mm
687G	C687N	Chromic Gut 5-0 18" P-3 Fine Reverse Cutting 3/8C 13mm
690G	A690N	Nylon 5-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
691G	A691N	Nylon 4-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
697H	A697N	Nylon 6-0 Black 18" P-1 Fine Reverse Cutting 3/8C 11mm
698H	A698N	Nylon 5-0 Black 18" P-3 Fine Reverse Cutting 3/8C 13mm
699G	A699N	Nylon 4-0 Black 18" P-3 Fine Reverse Cutting 3/8C 13mm
707G	D707N	Silk, 6-0, Black, 18", DR13, 3/8 Circle Taper Point 13mm
770G	B770N	Plain Gut 6-0 18" G-1 Fine Reverse Cutting 3/8C 11mm Double Armed
772G	B772N	Plain Gut 5-0 18" G-3 Fine Reverse Cutting 3/8C 13mm Double Armed
7733G	D7733N	Silk, 7-0, Black, 18", DSL7 DSL7, 3/8 Circle Side Cutting Lancet 6.6mm

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Long Description
774G	B774N	Plain Gut 6-0 18" G-1 Fine Reverse Cutting 3/8C 11mm
780G	D780N	Silk 6-0 Black 18" G-1 Fine Reverse Cutting 3/8C 11mm Double Armed
783G	D783N	Silk 4-0 Black 18" G-3 Fine Reverse Cutting 3/8C 13mm Double Armed
790G	C790N	Chromic Gut 6-0 18" G-1 Fine Reverse Cutting 3/8C 11mm Double Armed
792G	C792N	Chromic Gut 5-0 18" G-3 Fine Reverse Cutting 3/8C 13mm Double Armed
793G	C793N	Chromic Gut 4-0 18" G-3 Fine Reverse Cutting 3/8C 13mm Double Armed
796G	C796N	Chromic Gut 6-0 18" G-1 Fine Reverse Cutting 3/8C 11mm
8521H	J8521D	Polypropylene 4-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed
8522H	J8522D	Polypropylene 3-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed
8523H	J8523D	Polypropylene 2-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed
8526H	J8526D	Polypropylene 4-0 Blue 30" SH-1 Taperpoint 1/2C 22mm Double Armed
8533H	J8533D	Polypropylene 2-0 Blue 48" SH Taperpoint 1/2C 26mm Double Armed
8534H	J8534D	Polypropylene 3-0 Blue 48" SH Taperpoint 1/2C 26mm Double Armed
8554H	J8554D	Polypropylene 5-0 Blue 24" RB-2 Taperpoint 1/2C 13mm Double Armed
8555H	J8555D	Polypropylene 5-0 Blue 24" RB-1 Taperpoint 1/2C 17mm Double Armed
8556H	J8556D	Polypropylene 5-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8557H	J8557D	Polypropylene 4-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8558H	J8558D	Polypropylene 3-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8580H	J8580D	Polypropylene 5-0 Blue 36" BB Taperpoint 3/8C 17mm Double Armed
8581H	J8581D	Polypropylene 4-0 Blue 36" BB Taperpoint 3/8C 17mm Double Armed
8617G	J8617N	Polypropylene 6-0 Blue 18" PC-1 Fine Conventional Cutting 3/8C 13mm
8618G	J8618N	Polypropylene 5-0 Blue 18" PC-1 Fine Conventional Cutting 3/8C 13mm
8623H	J8623N	Polypropylene, 2-0, Blue, 30", GS60, Straight Straight Cutting 60mm
8631G	J8631N	Polypropylene 4-0 Blue 18" PC-5 Fine Conventional Cutting 3/8C 19mm
8632G	J8632N	Polypropylene 3-0 Blue 18" PC-5 Fine Conventional Cutting 3/8C 19mm
8634G	J8634N	Polypropylene 4-0 Blue 18" PC-3 Fine Conventional Cutting 3/8C 16mm
8635G	J8635N	Polypropylene 5-0 Blue 18" PC-3 Fine Conventional Cutting 3/8C 16mm
8661G	J8661N	Polypropylene 5-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8663G	J8663N	Polypropylene 3-0 Blue 18" PS-1 Fine Reverse Cutting 3/8C 24mm
8665G	J8665N	Polypropylene 3-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8680G	J8680N	Polypropylene 6-0 Blue 18" PS-3 Fine Reverse Cutting 3/8C 16mm
8681G	J8681N	Polypropylene 5-0 Blue 18" PS-3 Fine Reverse Cutting 3/8C 16mm
8682G	J8682N	Polypropylene 4-0 Blue 18" PS-2 Fine Reverse Cutting 3/8C 19mm
8683G	J8683N	Polypropylene 4-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8684G	J8684N	Polypropylene 3-0 Blue 18" FS-1 Reverse Cutting 3/8C 24mm
8685H	J8685N	Polypropylene 2-0 Blue 18" FS Reverse Cutting 3/8C 26mm
8686G	J8686N	Polypropylene 5-0 Blue 18" PS-2 Fine Reverse Cutting 3/8C 19mm

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Long Description
8687H	J8687N	Polypropylene 3-0 Blue 18" PS-2 Fine Reverse Cutting 3/8C 19mm
8695G	J8695N	Polypropylene 6-0 Blue 18" P-3 Fine Reverse Cutting 3/8C 13mm
8697G	J8697N	Polypropylene 6-0 Blue 18" P-1 Fine Reverse Cutting 3/8C 11mm
8698G	J8698N	Polypropylene 5-0 Blue 18" P-3 Fine Reverse Cutting 3/8C 13mm
8699G	J8699N	Polypropylene 4-0 Blue 18" P-3 Fine Reverse Cutting 3/8C 13mm
8701H	J8701D	Polypropylene 7-0 Blue 18" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8702H	J8702D	Polypropylene 7-0 Blue 24" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8703H	J8703D	Polypropylene 7-0 Blue 30" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8704H	J8704D	Polypropylene 7-0 Blue 24" CC Tapercut 3/8C 9.3mm Double Armed
8706H	J8706D	Polypropylene 6-0 Blue 30" C-1 Taperpoint 3/8C 13mm Double Armed
8707H	J8707D	Polypropylene 6-0 Blue 30" CC-1 Tapercut Cardio 3/8C 13mm Double Armed
8709H	J8709D	Polypropylene 6-0 Blue 30" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8710H	J8710D	Polypropylene 5-0 Blue 30" RB-2 Taperpoint 1/2C 13mm Double Armed
8711H	J8711D	Polypropylene 6-0 Blue 30" RB-2 Taperpoint 1/2C 13mm Double Armed
8713H	J8713D	Polypropylene 5-0 Blue 18" RB-2 Taperpoint 1/2C 13mm Double Armed
8714H	J8714D	Polypropylene 6-0 Blue 18" RB-2 Taperpoint 1/2C 13mm Double Armed
8717H	J8717D	Polypropylene 5-0 Blue 18" C-1 Taperpoint 3/8C 13mm Double Armed
8718H	J8718D	Polypropylene 6-0 Blue 18" C-1 Taperpoint 3/8C 13mm Double Armed
8720H	J8720D	Polypropylene 5-0 Blue 36" C-1 Taperpoint 3/8C 13mm Double Armed
8725H	J8725D	Polypropylene 5-0 Blue 24" C-1 Taperpoint 3/8C 13mm Double Armed
8726H	J8726D	Polypropylene 6-0 Blue 24" C-1 Taperpoint 3/8C 13mm Double Armed
8756H	J8756D	Polypropylene 5-0 Blue 18" RB-1 Taperpoint 1/2C 17mm Double Armed
8762H	J8762D	Polypropylene 3-0 Blue 30" SH-1 Taperpoint 1/2C 22mm Double Armed
8805H	J8805D	Polypropylene 6-0 Blue 24" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8806H	J8806D	Polypropylene 6-0 Blue 18" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8807H	J8807D	Polypropylene 6-0 Blue 24" CC Tapercut 3/8C 9.3mm Double Armed
8842H	J8842D	Polypropylene 3-0 Blue 36" MH Taperpoint 1/2C 36mm Double Armed
8843H	J8843D	Polypropylene 2-0 Blue 36" MH Taperpoint 1/2C 36mm Double Armed
8853H	J8853D	Polypropylene 2-0 Blue 48" MH Taperpoint 1/2C 36mm Double Armed
8935H	J8935D	Polypropylene 4-0 Blue 36" V-5 Tapercut Cardio 1/2C 17mm Double Armed
8975H	J8975D	Polypropylene 4-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
8976H	J8976D	Polypropylene 3-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
8977H	J8977D	Polypropylene 2-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
G182H	C182N	Chromic Gut, 3-0, 30", HR22, 1/2 Circle Taper Point 22mm
G667G	AG667N	Nylon 4-0 Blue 18" PS-2 Fine Reverse Cutting 3/8C 19mm
H821H	B821N	Plain Gut 4-0 27" FS-2 Reverse Cutting 3/8C 19mm
H822H	B822N	Plain Gut 3-0 27" FS-2 Reverse Cutting 3/8C 19mm

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Long Description
J232H	G232N	PolySyn™ , 3-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
J232H	G332N	PolySyn™ , 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
J232H	G416N	PolySyn™ , 3-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
J260H	G260N	PolySyn™ , 0, Undyed, 27", HR36, 1/2 Circle Taper Point 36mm
J266H	G266N	PolySyn 2-0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J267H	G267N	PolySyn 0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J268H	G268N	PolySyn 1 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J269H	G269N	PolySyn™ , 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
J310H	G310N	PolySyn™ , 4-0, Violet, 27", DGL19, 1/2 Circle Taper Point 20mm
J311H	G311N	PolySyn™ , 3-0, Violet, 27", DGL19, 1/2 Circle Taper Point 20mm
J316H	G316N	PolySyn™ , 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
J317H	G317N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
J333H	G333N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
J339H	G339N	PolySyn™ , 2-0, Violet, 27", HR36, 1/2 Circle Taper Point 36mm
J340H	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
J341H	G341N	PolySyn™ , 1, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
J345H	G354N	PolySyn™ , 2-0, Violet, 36", HRX36, 1/2 Circle Taper Point 36.6mm
J352H	G352N	PolySyn™ , 0, Violet, 27", HRX40, 1/2 Circle Taper Point 40mm
J358H	G358N	PolySyn™ , 0, Violet, 36", HRX40, 1/2 Circle Taper Point 40mm
J359H	G359N	PolySyn™ , 1, Violet, 36", HR40, 1/2 Circle Taper Point 40mm
J392H	G392N	PolySyn 4-0 Violet 18" FS-2 Reverse Cutting 3/8C 19mm
J393H	G393N	PolySyn 3-0 Violet 18" FS-2 Reverse Cutting 3/8C 19mm
J422H	G422N	PolySyn 4-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
J423H	G423N	PolySyn 3-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
J426H	G426N	PolySyn 4-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
J427H	G427N	PolySyn 3-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
J428H	G428N	PolySyn 2-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
J441H	G441N	PolySyn 4-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J442H	G442N	PolySyn 3-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J443H	G443N	PolySyn 2-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J463G	G463N	PolySyn 5-0 Violet 18" P-3 Fine Reverse Cutting 3/8C 13mm
J464G	G464N	PolySyn 4-0 Violet 18" P-3 Fine Reverse Cutting 3/8C 13mm
J489G	G489N	PolySyn 6-0 Undyed 18" P-1 Fine Reverse Cutting 3/8C 11mm
J492G	G492N	PolySyn 6-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
J493G	G493N	PolySyn 5-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
J494G	G494N	PolySyn 4-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Long Description
J495G	G495N	PolySyn 5-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
J496G	G496N	PolySyn 4-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
J497G	G497N	PolySyn 3-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
J500G	G500N	PolySyn 5-0 Undyed 18" PS-3 Fine Reverse Cutting 3/8C 16mm
J507G	G507N	PolySyn 4-0 Undyed 18" PS-4 Fine Reverse Cutting 1/2C 16mm
J535H	G535N	PolySyn 1 Undyed 27" OS-6 Reverse Cutting 1/2C 36mm
J544G	G544N	PolySyn™, 6-0, Violet, 18", VSL7 VSL7, 1/4 Circle Side Cutting Lancet 6.6mm
J556G	G556N	PolySyn™, 6-0, Violet, 12", VSPO8 VSPO8, 1/4 Circle Side Cutting Lancet 8mm
J670G	G670N	PolySyn™, 6-0, Undyed, 18", VCL8 VCL8, 1/4 Circle Side Cutting Lancet 8.5mm
J683G	G683N	PolySyn 3-0 Undyed 18" PS-1 Fine Reverse Cutting 3/8C 24mm
J824G	G824N	PolySyn 3-0 Undyed 18" PC-5 Fine Conventional Cutting 3/8C 19mm
J839D	G839N	PolySyn™, 2-0, Undyed, 18", HR36, 1/2 Circle Taper Point 36mm
J845G	G845N	PolySyn 4-0 Undyed 18" PC-3 Fine Conventional Cutting 3/8C 16mm
J869H	G869N	PolySyn 2-0 Undyed 27" CP-2 Reverse Cutting 1/2C 26mm
J917H	M917N	PolySyn FA™, 2-0, Undyed, 36", HRT36, 1/2 Circle Taper Cutting 36mm
J935H	G935N	PolySyn 4-0 Undyed 27" PS-1 Fine Reverse Cutting 3/8C 24mm
J936H	G936N	PolySyn 3-0 Undyed 27" PS-1 Fine Reverse Cutting 3/8C 24mm
J945H	G945N	PolySyn™, 2-0, Undyed, 36", HR36, 1/2 Circle Taper Point 36mm
K832H	D832N	Silk, 3-0, Black, 30", HR26, 1/2 Circle Taper Point 26mm
R633H	F633N	Polyviolene™, 4-0, White, 18", DS18, 3/8 Circle Reverse Cutting 19mm
R647H	F647D	Polyviolene™, 3-0, White, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
R665H	F665D	Polyviolene™, 2-0, White, 18", DS26, 3/8 Circle Reverse Cutting 26mm
VR214	M214N	PolySyn FA™, 4-0, Undyed, 27", HR17, 1/2 Circle Taper Point 18mm
VR426	M426N	PolySyn FA 4-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
VR493	M493N	PolySyn FA 5-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
VR494	M494N	PolySyn FA 4-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
VR496	M496N	PolySyn FA 4-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
VR497	M497N	PolySyn FA 3-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
X513H	EX513D	Polyviolene 3-0 Green 36" SH Taperpoint 1/2C 26mm Double Armed
X519H	E519N	Polyviolene 2 Green 30" OS-4 Reverse Cutting 1/2C 22mm
X523H	EX523D	Polyviolene 2-0 White 36" SH Taperpoint 1/2C 26mm Double Armed
X551H	EX551D	Polyviolene 4-0 Green 30" RB-1 Taperpoint 1/2C 17mm Double Armed
X556H	EX556D	Polyviolene 5-0 Green 36" RB-1 Taperpoint 1/2C 17mm Double Armed
X557H	EX557D	Polyviolene 4-0 Green 36" RB-1 Taperpoint 1/2c 17mm Double Armed
X558H	EX558D	Polyviolene 3-0 Green 36" RB-1 Taperpoint 1/2C 17mm Double Armed
X563H	EX563D	Polyviolene 2-0 Green 30" SH Taperpoint 1/2C 26mm Double Armed

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Long Description
X762H	EX762D	Polyviolene 3-0 Green 36" SH-1 Taperpoint 1/2C 22mm Double Armed
X763H	EX763D	Polyviolene 2-0 Green 36" SH-1 Taperpoint 1/2C 22mm Double Armed
X927H	EX927D	Polyviolene 2-0 Green 30" V-4 Taper Cutting 3/8C 17mm Double Armed
X977H	EX977D	Polyviolene 2-0 Green 30" V-7 Taper Cutting 1/2C 26mm Double Armed
Y266H	Y266N	Monoderm 2-0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
Y426H	Y426N	Monoderm 4-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
Y427H	Y427N	Monoderm 3-0 Undyed 27" PS-2 Fine Reverse Cutting 3/8C 19mm
Y432H	Y432N	Monoderm™, 6-0, Undyed, 27", HR13, 1/2 Circle Taper Point 13mm
Y433H	Y433N	Monoderm™, 5-0, Undyed, 27", HR13, 1/2 Circle Taper Point 13mm
Y434H	Y434N	Monoderm™, 4-0, Undyed, 27", HR13, 1/2 Circle Taper Point 13mm
Y463G	Y463N	Monoderm 5-0 Violet 18" P-3 Fine Reverse Cutting 3/8C 13mm
Y464G	Y464N	Monoderm 4-0 Violet 18" P-3 Fine Reverse Cutting 3/8C 13mm
Y492G	Y492N	Monoderm 6-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
Y493G	Y493N	Monoderm 5-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
Y494G	Y494N	Monoderm 4-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
Y495G	Y495N	Monoderm 5-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Y496G	Y496N	Monoderm 4-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Y497G	Y497N	Monoderm 3-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Y513G	Y513N	Monoderm 4-0 Violet 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Y682H	Y682N	Monoderm 4-0 Undyed 18" PS-1 Fine Reverse Cutting 3/8C 24mm
Y823G	Y823N	Monoderm 4-0 Undyed 18" PC-5 Fine Conventional Cutting 3/8C 19mm
Y834G	Y834N	Monoderm 5-0 Undyed 18" PC-1 Fine Conventional Cutting 3/8C 13mm
Y844G	Y844N	Monoderm 5-0 Undyed 18" PC-3 Fine Conventional Cutting 3/8C 16mm
Y845G	Y845N	Monoderm 4-0 Undyed 18" PC-3 Fine Conventional Cutting 3/8C 16mm
Y935H	Y935N	Monoderm 4-0 Undyed 27" PS-1 Fine Reverse Cutting 3/8C 24mm
Y936H	Y936N	Monoderm 3-0 Undyed 27" PS-1 Fine Reverse Cutting 3/8C 24mm
Z305H	R305D	PDO, 3-0, Violet, 27", 1/2 Circle Taper Point 17mm
Z422H	R422N	PDO 4-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
Z423H	R423N	PDO 3-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
Z443H	R443N	PDO 2-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
Z493G	R493N	PDO 5-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
Z494G	R494N	PDO 4-0 Undyed 18" P-3 Fine Reverse Cutting 3/8C 13mm
Z496G	R496N	PDO 4-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Z497G	R497N	PDO 3-0 Undyed 18" PS-2 Fine Reverse Cutting 3/8C 19mm
Z513G	R513N	PDO 4-0 Violet 18" PS-2 Fine Reverse Cutting 3/8C 19mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
3CL812	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
CC883	Y417D	Monoderm™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
CG812	C128N	Chromic Gut, 0, 30", HR36, 1/2 Circle Taper Point 36mm
CG883	C123N	Chromic Gut, 2-0, 30", HR26, 1/2 Circle Taper Point 26mm
CL2M	G260N	PolySyn™ , 0, Undyed, 27", HR36, 1/2 Circle Taper Point 36mm
CL461	G311N	PolySyn™ , 3-0, Violet, 27", DGL19, 1/2 Circle Taper Point 20mm
CL4M	G339N	PolySyn™ , 2-0, Violet, 27", HR36, 1/2 Circle Taper Point 36mm
CL5M	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
CL71M	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
CL802	G352N	PolySyn™ , 0, Violet, 27", HRX40, 1/2 Circle Taper Point 40mm
CL811	G339N	PolySyn™ , 2-0, Violet, 27", HR36, 1/2 Circle Taper Point 36mm
CL812	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
CL829	G232N	PolySyn™ , 3-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
CL829	G332N	PolySyn™ , 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
CL830	G269N	PolySyn™ , 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
CL844	G260N	PolySyn™ , 0, Undyed, 27", HR36, 1/2 Circle Taper Point 36mm
CL882	G316N	PolySyn™ , 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
CL883	G317N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
CL883	G333N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
CL924	G0518N	PolySyn™ , 0, Violet, 36", HRX36, 1/2 Circle Taper Point 36.6mm
CL930	G339N	PolySyn™ , 2-0, Violet, 27", HR36, 1/2 Circle Taper Point 36mm
CL931	G340N	PolySyn™ , 0, Violet, 27", HRX36, 1/2 Circle Taper Point 36.6mm
CL971	G0518N	PolySyn™ , 0, Violet, 36", HRX36, 1/2 Circle Taper Point 36.6mm
CM830	Y417D	Monoderm™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
CP844	J8444N	Polypropylene, 0, Blue, 30", HRT36, 1/2 Circle Taper Cutting 36mm
CS745	D2543N	Silk, 1, White, 30", HR36, 1/2 Circle Taper Point 36mm
GC123	Y417D	Monoderm™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
GG121	C121N	Chromic Gut, 4-0, 30", HR26, 1/2 Circle Taper Point 26mm
GG122	C122N	Chromic Gut, 3-0, 30", HR26, 1/2 Circle Taper Point 26mm
GG123	C123N	Chromic Gut, 2-0, 30", HR26, 1/2 Circle Taper Point 26mm
GG128	C128N	Chromic Gut, 0, 30", HR36, 1/2 Circle Taper Point 36mm
GL122	G316N	PolySyn™ , 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL123	G317N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL123	G333N	PolySyn™ , 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL127	G339N	PolySyn™ , 2-0, Violet, 27", HR36, 1/2 Circle Taper Point 36mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
GL322	G232N	PolySyn™, 3-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
GL322	G332N	PolySyn™, 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL323	G269N	PolySyn™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
GL66M	G317N	PolySyn™, 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL66M	G333N	PolySyn™, 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GL67M	G316N	PolySyn™, 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GLS122	G316N	PolySyn™, 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GLS123	G317N	PolySyn™, 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GLS123	G333N	PolySyn™, 2-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GLS322	G232N	PolySyn™, 3-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
GLS322	G332N	PolySyn™, 3-0, Violet, 27", HR26, 1/2 Circle Taper Point 26mm
GLS323	G269N	PolySyn™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
GM323	Y417D	Monoderm™, 2-0, Undyed, 27", HR26, 1/2 Circle Taper Point 26mm
SC5587G	Y463N	Monoderm™, 5-0, Violet, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SC5590G	Y464N	Monoderm™, 4-0, Violet, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SC5616G	Y495N	Monoderm™, 5-0, Undyed, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SC5617G	Y496N	Monoderm™, 4-0, Undyed, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SC5618G	Y497N	Monoderm™, 3-0, Undyed, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SC5627G	Y513N	Monoderm™, 4-0, Violet, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SC5637G	Y426N	Monoderm™, 4-0, Undyed, 27", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SC5638G	Y427N	Monoderm™, 3-0, Undyed, 27", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SC5679	Y682N	Monoderm™, 4-0, Undyed, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
SC5688G	Y492N	Monoderm™, 6-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SC5689G	Y493N	Monoderm™, 5-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SC5690G	Y494N	Monoderm™, 4-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SG1645	B770N	Plain Gut, 6-0, 18", DSM11 DSM11, 3/8 Circle Precision Reverse Cutting 11mm
SG5163G	C1638N	Chromic Gut, 3-0, 27", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SG5163G	C636N	Chromic Gut, 3-0, 30", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SG5637	C1637N	Chromic Gut, 4-0, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SG5637G	C1637N	Chromic Gut, 4-0, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SG5644	C1654N	Chromic Gut, 4-0, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SG5644G	C1654N	Chromic Gut, 4-0, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SG5654	B1644N	Plain Gut, 4-0, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SG5686	B686N	Plain Gut, 5-0, 18", DSM12, 3/8 Circle Precision Reverse Cutting 12mm
SM5624	Y495N	Monoderm™, 5-0, Undyed, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
SM5627	Y496N	Monoderm™, 4-0, Undyed, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SM5628	Y497N	Monoderm™, 3-0, Undyed, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SM5637	Y426N	Monoderm™, 4-0, Undyed, 27", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SM5638	Y427N	Monoderm™, 3-0, Undyed, 27", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SM5679	Y682N	Monoderm™, 4-0, Undyed, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
SM5687	Y493N	Monoderm™, 5-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SM5688	Y492N	Monoderm™, 6-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SM5690	Y494N	Monoderm™, 4-0, Undyed, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SN1964	A1864N	Nylon, 4-0, Black, 18", DSCM16, 3/8 Circle Precision Conventional Cutting 16mm
SN3686	A1696N	Nylon, 7-0, Black, 18", DSM11, 3/8 Circle Precision Reverse Cutting 11mm
SN5662	A1662N	Nylon, 4-0, Black, 18", HSM16, 1/2 Circle Precision Reverse Cutting 16mm
SN5662G	A1662N	Nylon, 4-0, Black, 18", HSM16, 1/2 Circle Precision Reverse Cutting 16mm
SN5663	A1663N	Nylon, 3-0, Black, 18", PS-1, 3/8 Circle Precision Reverse Cutting 24mm
SN5663G	A1663N	Nylon, 3-0, Black, 18", PS-1, 3/8 Circle Precision Reverse Cutting 24mm
SN5665G	A1281N	Nylon, 6-0, Black, 18", DSM18, 3/8 Circle Precision Reverse Cutting 18mm
SN5666	A1666N	Nylon, 5-0, Black, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SN5666G	A1666N	Nylon, 5-0, Black, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SN5667	A1667N	Nylon, 4-0, Black, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SN5667G	A1667N	Nylon, 4-0, Black, 18", PS-2, 3/8 Circle Precision Reverse Cutting 19mm
SN5668	A0556N	Nylon, 5-0, Black, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SN5668	A1668N	Nylon, 5-0, Black, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SN5668G	A0556N	Nylon, 5-0, Black, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SN5668G	A1668N	Nylon, 5-0, Black, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SN5669	A1669N	Nylon, 3-0, Black, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SN5669G	A1669N	Nylon, 3-0, Black, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SN5696	A1698N	Nylon, 6-0, Black, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SN5696G	A1698N	Nylon, 6-0, Black, 18", P-3, 3/8 Circle Precision Reverse Cutting 13mm
SN5862G	A1629N	Nylon, 4-0, Black, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
SP1634	J8634N	Polypropylene, 4-0, Blue, 18", DSCM16, 3/8 Circle Precision Conventional Cutting 16mm
SP1635	J8635N	Polypropylene, 5-0, Blue, 18", DSCM16, 3/8 Circle Precision Conventional Cutting 16mm
SP5663	J8663N	Polypropylene, 3-0, Blue, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
SP5663G	J8663N	Polypropylene, 3-0, Blue, 18", DSM24, 3/8 Circle Precision Reverse Cutting 24mm
SP5680G	J8680N	Polypropylene, 6-0, Blue, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SP5681	J8681N	Polypropylene, 5-0, Blue, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm
SP5681G	J8681N	Polypropylene, 5-0, Blue, 18", DSM16, 3/8 Circle Precision Reverse Cutting 16mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
SP5686	J8661N	Polypropylene, 5-0, Blue, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SP5686G	J8661N	Polypropylene, 5-0, Blue, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SP5687	J8665N	Polypropylene, 3-0, Blue, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SP5687G	J8665N	Polypropylene, 3-0, Blue, 18", DSM19, 3/8 Circle Precision Reverse Cutting 19mm
SP5695	J8617N	Polypropylene, 6-0, Blue, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SP5695G	J8617N	Polypropylene, 6-0, Blue, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SP623	J8623N	Polypropylene, 2-0, Blue, 30", GS60, Straight Straight Cutting 60mm
SP623BC	J8623N	Polypropylene, 2-0, Blue, 30", GS60, Straight Straight Cutting 60mm
SP629	J8629N	Polypropylene, 4-0, Blue, 18", DS24, 3/8 Circle Reverse Cutting 24mm
SS1639G	D639N	Silk, 6-0, Black, 18", DSM11, 3/8 Circle Precision Reverse Cutting 11mm
SS5639	D639N	Silk, 6-0, Black, 18", DSM11, 3/8 Circle Precision Reverse Cutting 11mm
SS5640	D640N	Silk, 5-0, Black, 18", DSM13, 3/8 Circle Precision Reverse Cutting 13mm
SS623	D623N	Silk, 2-0, Black, 30", GS60, Straight Straight Cutting 60mm
UM886	Y432N	Monoderm™, 6-0, Undyed, 27", HR13, 1/2 Circle Taper Point 13mm

N=box of 12; D=box of 36



Surgical Specialties Corporation, 247 Station Drive, Suite NE1, Westwood, MA 02090, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are property of their respective owners. ©2018 Surgical Specialties Corporation.
All Rights Reserved. SSC-028 R3 11/18

SURGICAL
SPECALTIES
CORPORATION™