





With full-wave rectified high frequency, the **HF100H+** is MinXray's most cost-effective units for medical imaging. The key feature for the **HF100H+** is the output of these units is independent of line voltage variations. Features that are included in the **HF100H+** include:

- CPU control of x-ray tube
- 5 memory stations for storing x-ray exposure techniques
- Automatic, dynamic line compensation
- Compact, lightweight
- Seamless integration with DR system is available



Where you care, we're there.

## HF100H+

CONSTRUCTION	l Mo	nobloc HF generator, full bridge inverter
оитрит	25 :	mA ® 40-60 kVDC, mA ® 62-80 kVDC, mA ® 82-100 kVDC,
KVDC STEPS	2K\	/ step
mAs	2.0	- 150
TIMER	Ulti 0.0	ra High Resolution 3-4.0 sec, 0.01 increments
X-RAY TUBE	Tosl	hiba D-124S
FOCAL SPOT	1.2	mm
ANODE HEAT ST	ORAGE 20	kHU
TOTAL FILTRATIO COLLIMATOR	ON INCL. 3.2	mm AI equivalent
INVERTER FREC	UENCY 601	(Hz
LINE VOLTAGE	AC :	Single Phase 100 - 140 V; 50/60 Hz
FUSE RATING	30 8	amp 125 kV
ELECTRICAL CONSUMPTION	20	amp at 125 kV
MAXIMUM CONSUMPTION		W @ 100 kVDC, 20 amp
SYSTEM NET WI	EIGHT 45.	5 lb (20kg)
TUBEHEAD DIM	IENSIONS 11.	5 x 9.1 x 17.7 in (29.2 x 23 x 45 cm)
TUBEHEAD WEI	GHT 44.	5 lb (20.2 kg)
COLLIMATOR	Con wit	ntinuously adjustable light beam type h central x-ray indicator
POWER CORD	20	ft (6.1 m)
EXPOSURE COR	D 8 ft	(2.44 m)



## **Ordering Information**

## HF100H+ High Frequency Portable X-Ray Unit

Variable kVDC (40-100), variable mA (30,25,20), complete with ultra high-resolution timer, adjustable light beam collimator, 5 memory stations, power cord, exposure switch and cord and technique chart.

## **XGS MKIV LW**

Gas Spring Portable Mobile Stand with 16 in wheels for **HF100H+** 

For more information about MinXray products, contact your local dealer or MinXray at:

MinXray, Inc. 3611 Commercial Avenue Northbrook, Illinois 60062 USA 1-847-564-0323 Toll Free: 1-800-221-2245 (U.S.&CA.)
E-mail: info@minxray.com
www.minxray.com



